



Splash! VPS Impression Material

Splash! is the impression material for the clinician who demands reliability and value in a VPS impression material. Splash! features excellent flow properties, exceptional tear strength and a wild berry scent that patients will love.

Available in two set times, five high-contrast viscosities and four delivery methods, Splash! is the superhydrophilic impression material for every technique and every impression.

- Super hydrophilic
- Two set times (2:15 and 4:30)
- Available in five high-contrast colors and four delivery methods
- Wild berry scent



Features	Benefits
Low contact angle	Increased hydrophilicity means more accurate castings
High tear strength	Will not tear upon removal from gingival sulcus
High dimensional stability	Multiple pours and exceptional detail reproduction indefinitely
Wild berry scent	Reduced gag reflex

Product	Item #
Splash! 2-Pack Cartridges Half-Time Set (2:15)	
Lite Body	SPD1231
Heavy Body	SPD1233
Contents: 2 - 48 ml cartridges, tips, instructions	

Product	Item #
Splash! Putty Paks	
Half-Time Set Putty Paks (2:45)	SPD1237
Contents: 12 - 18 ml paks, instructions	

Product	Item #
Splash! Putty Jars	
Regular Set Putty Jars (5:30)	SPD1212
Half-Time Set Putty Jars (2:45)	SPD1235
Contents: 1 - 250 ml jar base, 1 - 250 ml jar catalyst, mixing spoons, instructions	