



MEDESY[®]
ITALY



INSTRUCTIONS FOR USE — GAMMAFIX (EN)



EN New **GAMMAFIX** trays designed to facilitate, speed up and make the medical device reprocessing process safer. Their special shape allows surgical instruments to be handled in maximum safety during washing and sterilization procedures, thus considerably reducing the risk of accidental injuries for the user.

The trays are supplied with silicone rails to be secured; here are the simple illustrative images. The silicone rails can be secured in different positions, in accordance with the anatomy of the medical devices which will then be contained in the tray: it is recommended to find the ideal position so that the medical devices are held firmly and safely.

If necessary, cut the excess portion of the silicone rails with scissors.

