

Essential Expo Offers

Prices valid until
29th May 2023.

Meet with experts from your favorite brands
and the team from HealthCare Essentials.



Find Us
On The Red
Carpet!

- ☎ 0508 802 557
- 🌐 www.healthcareessentials.nz
- ✉ sales@healthcareessentials.nz
- 📷 [healthcareessentialsnz](https://www.instagram.com/healthcareessentialsnz)
- 📌 [HealthcareEssentialsDental](https://www.facebook.com/HealthcareEssentialsDental)

Prices exclude GST and freight. Prices may be subject to change.
Some items may be on backorder.



Auckland/North Shore/Northland	Emma-Jane	021 818 919
Auckland/Waitakere/Pukekohe	Fiona	021 922 014
Waikato/Bay of Plenty/Taupo	Vicki	021 899 729
Wellington/Lower NI/Gisborne	Anne	021 450 655
South Island	Rob	021 796 360

The Entire Range!
15% OFF

∞ BIOCLEAR



BT
MATRIX

* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.

centrix®

MAKING DENTISTRY EASIER®



FluoroDose®

5% Sodium Fluoride Varnish with Xylitol



DENTAL ADVISOR TOP AWARD:
FLUORIDE VARNISH
6 YEARS RUNNING.



Each unit contains:
one single dose (0.3ml) of FluoroDose
Plus one clean Benda® Brush applicator.

Essential Expo Offer

120 & 1200 Packs

Buy 4
Get 1 FREE*

120 Pack \$320.00
1200 Pack \$2,050.00



Mixed 40 Pack

Buy 4
Get 1 FREE*

Mixed Pack \$110.00
8 of each flavour

ENCORE® D/C Mini Kits

Dual-Cure Core Buildup Composite

Kits include 12x single-dose .5mL cartridges & 24 Mixing Nozzles



Natural or
Contrast:

\$320.00

ACCESS® EDGE

High-Viscosity
Kaolin Clay-Based
Retraction
Paste



Contains 15 prefilled tubes,
30 dispensing tips, 10 large
GingiCaps and 10 regular
GingiCaps

\$284.00



TEMPIT® FAMILY

Unique Temporary Filling Materials.
One Delivery System.

30 x .35g
prefilled tips

Essential Self Cure Light Cure

\$180.00 \$175.00 \$210.00

Essential Expo Offer

\$199
Value

Buy 1 Get a
Snub Nose
Applicator
Gun FREE*



Essential Expo Offer

Buy 4 Get 1 FREE*

MAKING DENTISTRY EASIER®

* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.

ESSENTIAL EXPO OFFERS

3




Clinician's Choice

Designed for better dentistry.

A.S.A.P

All Surface Access Polishers



	Code	Usually	15% OFF
Indirect Starter Kit	CC-130918	\$179.00	\$152.15
Regular Starter kit	CC-130917	\$326.00	\$277.10
New Bulk Pack small pre polisher	CC-130927	\$395.00	\$335.75
New Bulk Pack large pre polisher	CC-130928	\$395.00	\$335.75
New Bulk Pack small final polisher	CC-130925	\$395.00	\$335.75
New Bulk Pack large final polisher	CC-130926	\$395.00	\$335.75

DualForce

Sectional Matrix System



	Code	Usually	15% OFF
Starter Kit	CC-209021	\$604.00	\$513.40

Template

Ultra Quick Matrix Material



	Code	Usually	15% OFF
4 Pack	CC-028904	\$215.00	\$182.75
12x50ml Bulk Pack	CC-018912	\$499.00	\$424.15

Evanesce

Nano-Enhanced Universal Restorative Composite



Cling2

Resin Optimized Temporary Cement



	Code	Usually	15% OFF
5ml	CC-158900	\$110.00	\$93.50
10ml	CC-158910	\$163.00	\$138.55

Affinity

Hydroactive Impression Material



* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.

DIRECTA

The Entire Range!
20% OFF

**Luxator®
Periotomes**




**physics®
forceps**



FENDERWEDGE®



PRACTIPAL® TRAYS



CONTACTEZ®



The Entire Range!
15% OFF

Virofex®



parkell®



* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.

The Entire Range!
20% OFF



* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.

Essential Expo Offers

NEW USER - 20% OFF
Bone Regeneration Material
 MAXIMUM OF 3

FREE BLUE M 3PK
OF ORAL GEL
 When you order \$1,000.00 of Ethoss



ethoss® is the result of many years clinical experience with Guided Bone Regeneration in dental surgery. The human body is amazing and wants to heal. ethoss® has been developed to work with it by using fully biocompatible materials to help true host bone regeneration. We prefer to say that it is a material that will quickly and easily turn into bone to help speed along dental procedures.

ethoss® provides a stable scaffold that bio-absorbs. It is completely free from any animal or human content.

It provides an environment rich in phosphate and calcium ions that are absorbed at the required rate to be replaced by new host bone. This true regeneration leaves only host bone. The patient's diseased or damaged tissue is returned to its original state for optimal support of both dentition and dental implants.

ethoss® is:

- Synthetic bi-phasic paste with a simple syringe delivery system
- Sets in vivo to provide a highly stable graft
- No need to use an additional collagen membrane
- Fast conversion to bone – human biopsies show 50% bone at 12 weeks
- Fully replaced by the patient's own bone
- Does not suffer wash-out in bloody sites
- Maintains long term bone volume


ethoss® can be used in the following clinical situations:

- Guided Bone Regeneration
- Periodontitis
- Buccal defects
- Sinus grafts
- Socket filling


EthOss Synthetic Bone Regeneration (3 x 0.25cc) **ET-T020230S** **\$500.00 ea**

EthOss Synthetic Bone Regeneration (3 x 0.5cc) **ET-T050530S** **\$580.00 ea**

EthOss Synthetic Bone Regeneration (3 x 1.0cc) **ET-T101030S** **\$695.00 ea**



EthOss Case Studies Facebook Group



* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.

The Entire Range!
20% OFF

GERMIPHENE®



The Entire Range!
15% OFF

Astek
 INNOVATIONS




The Entire Range!
10% OFF

Microbrush International®



* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.



LEVANTE was created and sculpted in the same way the artist shapes its work of art, finely removing the superfluous elements in search for the essence of pure form. The handle perfectly reflects **PININFARINA** passion as well as strong expertise in combining elegance of forms with technology, where «**innovation is indeed humanized by the design originality**»

LEVANTE[™]
by pininfarina



Maestrale Elevators

With Maestrale MEDESYS redefines the standards of a classic instrument, elevating it in every aspect.



* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.

The Entire Range!
20% OFF



Essential Expo Offer

Buy 5 Get 5 FREE*

- Beautiful Injectable X
- Beautiful II
- Beautiful XSL
- Beautiful II LS



BeautiBond Xtreme

NEW!
Product



THE TRUE UNIVERSAL ADHESIVE
 Etch, prime and bond in one step!



- For a wide range of applications and adherends with one bottle
- Bonds equally well to wet and dry dentine
- No filler and no HEMA - protection against hydrolysis
- Compatible with self-etch, total-etch and selective-etch techniques
- Simplified procedure for any adherend surface
- No additional primer needed



* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.



We care for healthy smiles



TePe Kids




TePe Compact Tuft



TePe Universal Care



TePe Interspace



ISO Size	Wire Size	 x6 Blister Pack	 x8 Soft Bag	 x25 Soft Bag
0	0.4mm	TP-137610	TP-138210	TP-132410
1	0.45mm	TP-137620	TP-138220	TP-132420
2	0.5mm	TP-137630	TP-138230	TP-132430
3	0.6mm	TP-137640	TP-138240	TP-132440
4	0.7mm	TP-117650	TP-118250	TP-112450
5	0.8mm	TP-117660	TP-118260	TP-112460
6	1.1mm	TP-117670	TP-118270	TP-112470
7	1.3mm	TP-117680	TP-118280	TP-112480
8	1.5mm	-	TP-118290	-
Mixed	-	-	TP-112299	-



A choice that matters

By replacing fossil raw material with fossil-free pine oils, we have reduced the carbon footprint of this interdental brush by 80% – without compromising on product quality, efficiency or design.



* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.

The Entire Range!
15% OFF

VA
 VISTA | APEX



Seamfree

Composite Wetting Resin & Lubricant

New!

Quick-Stat™ Gel

Our Quick-Stat™ Gel is a highly effective 20% ferric sulfate gel used for immediate and profound control of bleeding. This gel works on both hard and soft tissue and is highly effective for the control of bleeding during impressions, restorations, and crown and bridge procedures. Quick-Stat™ Gel will not flow onto surrounding tissue and quickly controls bleeding. Delivered with Stat-Flo™ applicator tips, this hemostatic agent is easy to use and incredibly effective.



Quick-Stat™ FREE

Quick-Stat™ FREE is the first clear hemostatic agent with shortened hemostasis time. Metal ions in hemostatic agents cause the precipitation of proteins. The positively charged metal ions bind with the negative charge of the blood/serum proteins and you get coagulation because the denatured proteins help accumulate platelets and fibrinogen to speed up clot formation.

Not only that, the proprietary surfactant in Quick-Stat™ FREE allows better penetration to more effectively plug cut capillaries. This is the Astringent Hybrid Technology that makes Quick-Stat™ FREE 20% faster than other clear hemostatic agents. For controlling bleeding, especially in the aesthetic zone, Quick-Stat™ FREE is the clear choice.



ENDO|ULTRA®

EndoUltra® Cordless Ultrasonic Activator: Best-In-Class Ultrasonic Activation of Irrigants

V-407702

Usually \$1,680.00 **\$1,428.00**



* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.

The Entire Range!
15% OFF



NeoPUTTY

Bioactive bioceramic premixed root and pulp treatment



NeoMTA[®] 2

Powder/Liquid root & pulp treatment material



NeoSEALER[™] Flo


Bioactive bioceramic sealer

R-Motion[®]

XP ENDO[®] solutions

RACE[®] EVO

The Entire Range!
15% OFF



FKG
swiss endo

NEW!
Endodontic File Range

* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.

The Entire Range!
10% OFF

inibsa

ARDANEST
with adrenaline
(epinephrine) 1:200,000
1.8 ml

Each cartridge (1.8 ml) contains 72 mg of articaine hydrochloride and 0.0072 mg of adrenaline (epinephrine, as tartrate).
Each ml contains 40 mg of articaine hydrochloride and 0.04 mg of adrenaline (epinephrine, as tartrate).

50 CARTRIDGES OF 1.8 mL SOLUTION FOR INJECTION / OROMUCOSAL USE

Ardanest 4% 1:200,000
Articaine HCl 4% + Epinephrine (tartrate)

IN-4670

Usually \$57.00

\$51.30

50 PACK of 1.8ml Cartridges

50% OFF

inibsa

ARDANEST
with adrenaline
(epinephrine) 1:100,000
1.8 ml

Each cartridge (1.8 ml) contains 72 mg of articaine hydrochloride and 0.0144 mg of adrenaline (epinephrine, as tartrate).
Each ml contains 40 mg of articaine hydrochloride and 0.08 mg of adrenaline (epinephrine, as tartrate).

50 CARTRIDGES OF 1.8 mL SOLUTION FOR INJECTION / OROMUCOSAL USE

Ardanest 4% 1:100,000
Articaine HCl 4% + Epinephrine (tartrate)

IN-4669

50 PACK of 1.8ml Cartridges

Expiry 31/10/23



\$28.50

20% OFF

SWEEP AG200
ULTRASONIC CLEANER
With Basket

LR-AG200

Usually \$2,295.00

\$1,836.00



The Entire Range!
10% OFF

Valisafe

Class 6 Self-Adhesive Indicators

MS-V3504065

Usually \$65.00

\$58.50



En867-5 Helix Kit

MS-V3502034

Usually \$257.00

\$231.30

* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.

optident LTD

WHITE DENTAL BEAUTY

NOVON

WHITE DENTAL BEAUTY
PROFESSIONAL DISSOLVING WHITENING STRIPS
6% HP MINT FLAVOUR

The Entire Range!
15% OFF

New!

Oxyfresh® Toothpaste & Mouthwash

The Entire Range!
20% OFF


PDT

Paradise Dental Technologies

The Entire Range!
35% OFF

* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.

The Entire Range!
15% OFF



PULPDENT®



The Entire Range!
10% OFF



RØNVIG



Scaling Tips

The Scorpion insert tip is coated with a 2 µm thick and strong titanium nitride layer.

To maintain efficiency and for patient comfort the insert has to be changed when the titanium coating has disappeared.



Essential Expo Offer
Buy 3 Tips Get 1 FREE*



* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.




Optima E+ BAE380R

Cordless Endodontic Unit with Apex Locator



◀ Simple control system for ease of use



▲ Easy to read display



◀ 360° rotatable head

A brushless endodontic motor system with inbuilt apex locator, reciprocating features and presets for many file types.

Inbuilt Apex Locator

The most notable feature of the BAE380R is its inbuilt apex locator. This tool provides integrated length determination, real-time display of file location in the canal, and automatic stop or reverse functions for when the file reaches the apical stop. The unit can also work as a standalone apex locator using the EAL setting.

Safety Glide Path (reciprocation) Mode

Using the BAE380R's SGP (Safety Glide Path) or reciprocation mode, the angles can be adjusted, whereby the forward angle must be equal to the reverse angle. Reciprocation allows the dentist greater flexibility in endodontic procedures and has been used extensively in endodontics for many years.

Brushless Motor

The brushless motor of the BAE380R not only makes the motor more cost effective, but it is also one of the quietest motors currently on the market, ensuring improved patient comfort with its low noise levels.

Optima E+
BAE380R

M2
CW 450rpm
2.0Ncm

NEW PRODUCT

\$1,360.00

Usually
~~\$1,600.00~~

BA-BAE380R

Specifications & Operation Modes:

Torque Range	0.4Ncm-5.0Ncm (4mNm ~ 50mNm)
Speed Range	100rpm~2500rpm
Wireless charging	
Frequency Range	112-205KHz
Operation Modes	EAL, CW, CCW, SGP, ATR
Customizable	
Memory programs	10
Inbuilt endo system presets	32
Battery specification	Lithium 3.7V /2000mAh
Input:	~100V-240V 50Hz/60Hz 0.4A Max
Output:	DC5V/1A
EAL	(Electronic Apex Locator)
CW	(Clockwise)
CCW	(Counter-Clockwise)
SGP	(Safety Glide Path)
ATR	(Adaptive Torque Reverse)



* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.



TRADESHOW SPECIALS

2023
NEW ZEALAND



2+1*
BUY 2 GET 1 FREE

9+5*
BUY 9 GET 5 FREE

BUY 1 BOX OF RIVA
RECEIVE A BONUS
20% DISCOUNT

BUY 8 BOXES OF RIVA
RECEIVE A BONUS
ULTRAMAT S

You can mix and match the
Glass Ionomer Capsules*



ionglass™
BIOMIMETIC FORMULA



CALL FREE 0508 802 557

Fax: 04 3873874
sales@healthcareessentials.nz
www.healthcareessentials.nz

YOUR LOCAL REPRESENTATIVE:

Auckland / North Shore / Northland
Emma-Jane Somerfield 021 818 919

Auckland / Waitakere / Pukekohe
Fiona Zuccani-Heard 021 922 014

Waikato / Bay of Plenty / Taupo
Vicki Slow 021 899 729

Wellington / Lower North Island / Gisborne
Anne Carter 021 450 655

South Island
Rob Evans 021 796 360

*All prices exclude GST. Offers are valid at exhibitions and seminars attended by SDI in 2023. Offers cannot be combined. All prices exclude GST. All claims for free goods must be made within 30 days from the end of promotion date. Free goods will be of equal or lesser value. The above-mentioned offers are not applicable on corporate, government accounts or those customers who have received discounted pricing. Offers are valid in New Zealand only.

* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.



ESSENTIAL EXPO OFFERS

SDI | **RIVA SELF CURE**
BULK FILL GLASS HYBRID RESTORATIVE SYSTEM



THE IONGLASS LONG TERM BREAKTHROUGH

When dentine is missing, use Riva Self Cure to replace it. It is the best dental material available today that virtually mimics dentine. No adhesive is required, and sensitivity is non-existent.

- Contains SDI's proprietary ionglass™, a high ion releasing, radiopaque, bioactive glass that remineralises and prevents secondary caries formation
- Very high fluoride releasing and recharging ability
- Choice of viscosity - High Viscosity formulation is flowable enough to inject and highly sculptable and packable
- Exceptional Handling; Non-sticky, non-slumping and sculptable formulation
- Easy Activation Capsule - 70% easier to activate
- Ideal materials for dentine replacement
- Radiopaque
- BPA & HEMA free



BOX OF 50 CAPSULES \$230+ GST

RIVA PROTECT

SELF-CURING GLASS IONOMER USED AS A SEALANT, LINER AND SURFACE PROTECTOR

- ACP nanotechnology - supercharges remineralisation
- Very high fluoride releasing and recharging ability to protect teeth from acid attacks
- Low viscosity for easy placement and pit/fissure penetration
- BPA & HEMA free
- Contains SDI's proprietary ionglass™, a high ion releasing, radiopaque, bioactive glass that remineralises and fights against secondary caries formation
- Works in a moist environment - indicated for partially erupted teeth and challenging patients
- Available in white and pink shades - pink shade is ideal for high visibility requirements
- Bioactive proprietary ionglass™ technology - high ion releasing

BOX OF 50 CAPSULES \$233 + GST

SDI | **RIVA LIGHT CURE**
LIGHT CURED RESIN REINFORCED
GLASS IONOMER RESTORATIVE MATERIAL



THE MULTI AWARDED EASY HANDLING GLASS IONOMER

Riva Light Cure has excellent resistance to common stains such as black tea and red wine.

- Contains SDI's proprietary ionglass™, a high ion releasing, radiopaque, bioactive glass that remineralises and fights against secondary caries formation
- Exceptional handling; easy to sculpt and non-sticky - available in both regular and high viscosity
- Unique, patented, orange light safe tip, prevents curing in the tip and significantly extends working time
- Fabulous aesthetics
- Easy Activation Capsule - 70% easier to activate
- Highest colour stability to tea & coffee
- Excellent marginal integrity
- Long term wear resistance due to high surface hardness
- High flexural strengths
- BPA free



BOX OF 50 CAPSULES \$287 + GST



RIVA SILVER

MERCURY FREE SELF-CURING
GLASS IONOMER RESTORATIVE MATERIAL

- Silver reinforced
- Ideal for abfraction lesion restorations
- Ideal for bulk filling
- Great for core build ups
- Easy handling
- No shrinkage
- High fluoride releasing and recharging ability
- Contains SDI's proprietary ionglass™, a high ion releasing, radiopaque, bioactive glass that remineralises and fights against secondary caries formation
- High wear resistance
- Bioactive proprietary ionglass™ technology - high ion releasing to assist remineralisation
- Very high radiopacity
- HEMA & BPA free

BOX OF 50 CAPSULES \$263 + GST



RIVA LIGHT CURE HV

* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.

The Entire Range!
15% OFF


pinkwave™

VA
VISTA | APEX

PINK IS THE NEW CURE



PinkWave™ was distinctly designed to enhance the efficacy of curing procedures. PinkWave™ is equipped with Vista Apex's patented QuadWave™ Technology which leverages four different wavelengths to ensure peak performance. Along with the standard blue light, PinkWave™ also employs red, near infrared (NIR) and UV light, which make up the only pink curing light on the market. This innovative approach to curing has never been achieved before and results in both doctor and patient benefits.

- Built in Transilluminator
- 3 curing modes — Boost, Standard and Ramp
- 1515 mW/cm² on Standard mode and 1720 mW/cm² on Boost mode

PinkWave™ is the next generation of curing lights, equipped with patented QuadWave™ Technology, offering decreased composite shrinkage and increased depth of cure. With the largest curing light area on the market, PinkWave™ ensures maximum coverage. QuadWave™ Technology leverages four wavelengths into one collimated pink beam: UV, Blue, Red, and Near-infrared.

- 12.1 mm Curing Area
- Decreased Shrinkage
- Up to 8.5mm depth of cure

CURING LIGHT KIT

Includes:

- 1 x Cordless Curing Light
- 1 x Induction Charging Base
- 1 x Power Cord
- 5 x Autoclavable Light Shields
- 100 x Barrier Sleeves
- 3 x Pair Protective Eyewear

V-90734

Usually \$3,000.00 **\$2,550.00**

* Entire Range Discount does not apply to Buy-Get offers - standard pricing applies.



ESSENTIAL EXPO OFFERS