

# Practice

Essentials **Dec** 2023 **Jan** 2024



- ☎ 0508 802 557
- 🌐 [www.healthcareessentials.nz](http://www.healthcareessentials.nz)
- ✉ [sales@healthcareessentials.nz](mailto:sales@healthcareessentials.nz)
- 📷 [healthcareessentialsnz](https://www.instagram.com/healthcareessentialsnz)
- 📘 [HealthcareEssentialsDental](https://www.facebook.com/HealthcareEssentialsDental)

Prices exclude GST and freight. Prices may be subject to change.  
Some items may be on backorder.



## HealthCare Essentials

Auckland/North Shore/Northland	<b>Emma-Jane</b>	021 818 919
Auckland/Waitakere/Pukekohe	<b>Fiona</b>	021 922 014
Waikato/Bay of Plenty/Taupo	<b>Vicki</b>	021 899 729
Wellington/Lower NI/Gisborne	<b>Anne</b>	021 450 655
South Island	<b>Nicky</b>	021 796 360

# Practice Essentials

Dec 2023 Jan 2024

2	Topical Anaesthetic
3	Syringes
4	Anaesthetic
5-6	Diamond Burs
7	Burs
8	Finishing & Polishing
9-10	Bur Stands
11-12	Matrix Bands
13-14	Wedges
15-16	Matrix Kits
17	Rubber Dam
18	Caries Indicators & Desensitiser
19	Etch
20	Adhesive & Bond
21-24	Composite & Glass Ionomer
25	Curing Light
26	Articulating Paper & Interproximal Strips
27-28	Finishers & Polishers
29	Fluoride Varnish
30-32	Cassettes
33-34	Forceps
35-36	Impression Material
37	Temporary Cements & Materials
38	Impression Trays, Adhesive & Alginate
39-40	Impression Trays
41-43	Tips
44	Prophy Paste
45	Instruments
46	Ultrasonic Scaling
47-48	Whitening
49-50	Toothbrushes
51-53	Infection Control
54	Cotton Rolls & Micro Applicators
55	Disposables
56	Sleeves & Barriers
57	Gloves, Masks & Luer Lock Syringes
58	Ultrasonic Cleaning
59	Sterilisation Pouches & Christmas Closure Dates

# ZAP Topical Anaesthetic Gel



Temporarily eliminating pain makes everyone's job easier including your patient's.

**Effective**

Minimises the pain in oral mucosal tissue arising from needle punctures, deep scaling procedures, prosthetic adjustments, clamp or crown placement, removal of primary teeth and suture removal. The synergy of benzocaine and tetracaine results in a fast and long lasting anaesthetic.

**No-gagging**

Zap can be used to reduce the pharyngeal (gag) reflex associated with placing dental materials like impression trays and x-ray film into the oral cavity.

**Convenient**

Award winning one touch jar eliminates the chances of cross-contamination and dispenses the right amount of gel each time, reducing waste. The frosted dispenser shows you how much product is left in the jar so you can reorder efficiently.

**Fast-acting**

Topical anaesthesia within 30 seconds. Topical anesthesia for a minimum of 10 –15 minutes.\*

**Taste**

Cool Mint flavour eliminates any bitter aftertaste.



Award Winning **One Touch** 50ml Applicator Jar



**GLUTEN FREE**

**Effective in 30 seconds** for up to 10-15 minutes\*

**Contains 18% benzocaine and 2% tetracaine**

**Zap is ideal for:**

**General**

- Pre-injection anaesthetic
- Gag reflex from impressions, X-rays
- Alleviate apprehension
- Retraction cord placement
- Clamp placement

**Endodontics**

- Pre-injection anaesthetic
- Lubricant with rubber dams
- On file or reamer as a lubricant and anaesthetic to help eliminate the pain of intrapulpal injection

**Periodontics**

- Pre-injection anaesthetic
- Localised gingival curettage
- Deep scaling
- Suture removal
- Post-op discomfort

**Oral Surgery**

- Pre-injection anaesthetic
- Dry socket
- Suture removal

**Orthodontics**

- Minor mouth irritation
- Banding a molar
- Minimise gag reflex for more accurate impressions
- Arch bar removal

**Paedodontics**

- Pre-injection anaesthetic
- Removal of mobile deciduous teeth

**Prosthodontics**

- Denture discomfort
- Minimise gag reflex for more accurate impressions



**\$67.50**

Usually \$75.00  
50ml Applicator Jar

GP-HZAP

50 g



**ZAP**  
Topical Anaesthetic Gel  
Cool Mint

For topical use only  
DIRECTIONS  
1. Press down on the applicator jar  
2. Apply the gel to the area to be anesthetized.  
CONTENTS: Benzocaine, Tetracaine  
KEEP OUT OF REACH OF CHILDREN  
GERMIPHENE  
Grand Rapids, MI  
1.800.255.9912

\*Data on file. | Please refer to product label for complete usage and safety information.

### ASPIJECT® Passive

Injection syringe for infiltration and block anaesthesia

- Made from high-grade stainless steel and heat-resistant engineering polymer
- Automatic (passive) aspiration occurs when a slight pressure on the plunger is eased
- Self-aspirating without pullback minimises the risk of needle tip movement
- Can be sterilised by any method
- Material and construction guaranteed for 5 years
- New superior surface treatment

RV-2055 Aspiject Passive 1.8ml (CR)

RV-2056 Aspiject Passive 1.8ml (TR)

RV-2060 Aspiject Passive 2.2ml (CR)

RV-2061 Aspiject Passive 2.2ml (TR)

**\$254.15**  
Usually \$299.00



**5 Year Warranty**

### ASPIJECT® Active

Injection syringe for infiltration and block anaesthesia

The syringe has an ergonomic design that provides optimal balance and a tactile control of the injection.

**\$255.00**

Usually \$300.00

RV-2075 1.8ml (TR)

RV-2080 2.2ml (TR)

Self-tapping spiral hook



**5 Year Warranty**

### ASPIJECT® SIS

Sterile Injection System

Specially designed for use during surgical procedures. Self-aspirating dental cartridge syringe for infiltration and block anaesthesia.

**\$232.00**

Usually \$273.00

RV-2090 1.8ml



3 pack 1.8ml Autoclavable Barrels

**\$51.80**

Usually \$61.00

RV-2028

**5 Year Warranty**

### ASPIJECT® Paroject

Injection syringe for intraligamental local anaesthesia.

PAROJECT® is an unthreatening penstyle intraligamental syringe with a controlled dosage (0.6ml) permitting a slow non-traumatic injection.

**\$589.00**

Usually \$693.00

RV-2018 1.8ml

RV-2030 2.2ml

Paroject parts are available on request.



**5 Year Warranty**

### ASPIJECT® Astra

Aspirating Syringe 2.2ml

Self aspirating syringe using a standard anaesthetic cartridge. Manufactured in stainless steel to the highest standards.

**\$186.15**

Usually \$219.00

RV-2065 2.2ml



Available while stocks last!

**5 Year Warranty**

# ARDANEST 4%



Articaine HCl 4% + Epinephrine (tartrate) 1:200,000 or 1:100,000

## Ardanest 4% Articaine HCl 4% + Epinephrine (tartrate)

### 1:200,000



1 mL of solution for injection contains articaine hydrochloride 40.00 mg and adrenaline (epinephrine, as tartrate) 0.005 mg.

One cartridge (1.8 mL) contains articaine hydrochloride 72.00 mg and adrenaline (epinephrine, as tartrate) 0.009 mg.

Each 1.8 mL cartridge contains: 0.71 mg of sodium.

### Duration of Action<sup>1</sup>:




**Pulpar Anaesthesia**  
**45-60 min.**



**Soft Tissue**  
**120-300 min.**

<sup>1</sup> Malamed, SF. - Handbook of local Anaesthesia. Fifth Edition. Elsevier Mosby, 2004. p.71-73. 130.



July 2024  
Expiry:

25% OFF  
**\$42.75**

Regular Dated Price  
**\$57.00**

IN-4670 4% 1:200,000 50 PACK of 1.8ml Cartridges

## Ardanest 4% Articaine HCl 4% + Epinephrine (tartrate)

### 1:100,000



1 mL of solution for injection contains articaine hydrochloride 40.00 mg and adrenaline (epinephrine, as tartrate) 0.01 mg.


One cartridge (1.8 mL) contains articaine hydrochloride 72.00 mg and adrenaline (epinephrine, as tartrate) 0.018 mg.

Each 1.8 mL cartridge contains: 0.71 mg of sodium.

### Duration of Action<sup>1</sup>:



**Pulpar Anaesthesia**  
**60-75 min.**



**Soft Tissue**  
**180-360 min.**

<sup>1</sup> Malamed, SF. - Handbook of local Anaesthesia. Fifth Edition. Elsevier Mosby, 2004. p.71-73. 130.



Regular Dated Price  
**\$57.00**

IN-4669 4% 1:100,000 50 PACK of 1.8ml Cartridges

### Articaine anaesthesia:

- Articaine has faster onset of action (~2 min.) and longer duration than lidocaine and mepivacaine. This is due to its:
  - high liposolubility, which allows it to move better through soft and hard tissues
  - high protein binding, which makes it more firmly bound in the nerve, resulting in longer-lasting blockade or anaesthesia
- Articaine is 1.5 times more potent than lidocaine
- Articaine offers intermediate duration of action
- Articaine is a potent vasodilator, equal to lidocaine. Both cause more vasodilation than mepivacaine and prilocaine
- Local anaesthetics with articaine are formulated in combination with a vasoconstrictor such as epinephrine (adrenaline). Epinephrine increases the duration of action but can have short-term cardiovascular effects
- Indications for both A100 and A200:
  - local, infiltrative, or conductive anaesthesia in simple and complex dental procedures (medium to long duration routine and surgical procedures)
- The anaesthetic efficacy of 4% articaine with 1:200.000 epinephrine is comparable to 4% articaine with 1:100.000



Diamond Points FG  
Regular Grit

Single \$8.70 10 Pack 40% OFF 100 Pack\* 50% OFF



Specially designed for cavity cutting and tooth preparation.  
Available in: 88 different shapes, some shapes are also in short shank (categorised under 1000 class).

Item Code (SHDP+)	401M	401	301	301M	201S	201	101S	101	101L	202	202S	103	301CR	202CR	201CR	101CR	103CR	206CR	207CR	201SR	201R	104R	101R
Short	1401	1301			1201		1101	1101L	1202	3202	1103							0825	0826		1201R	1104R	1101R
SuperShort	3401	3301																			3201R	3104R	
Robot Points																							
ISO#	170/010	170/011	170/016	170/013	171/014	171/016	173/014	173/016	173/018	170/021	170/017	172/022	544/018	544/021	545/018	553/019	553/023	544/027	544/031	222/014	222/016	193/016	199/016
L(mm)	3	4	4.5	4	7	7	10	10	10	5	5	9	4	5	7	10	10	3.5	3.5	7	7	10	10
Item Code (SHDP+)	101SR	102R	103R	105R	106RD	107RD	201Ra	104Ra	101Ra	304	204	104	104L	411	411M	311	211	214	311RM	212R	213R	114S	114
Short		1102R	1103R								1204	1104		1411		1311	1211						
SuperShort		3102R			0836	0839					3204	3104											
Robot Points																							
ISO#	199/013	223/017	223/022	223/022	198/021	198/023	218/016	220/016	220/016	164/016	165/013	166/016	166/018	109/011	109/010	109/013	110/014	156/012	139/013	288/010	289/012	158/010	158/012
L(mm)	9	9	9	11.5	9	9	7	10	10	6	8	10	10	4	3	4.5	7	5	4.5	6	8	8	8
Item Code (SHDP+)	316	215S	215	216	115	116	421M	421	321M	321	120	421SS	421S	321S	435	335	150	440SS	440S	440	340	340S	461R
Short			1215	1216								1421SS	1421S					3440SS		1440	1340		3461R
SuperShort			3215									3421SS					0899			3440			
Robot Points																							
ISO#	247/013	249/012	249/014	249/017	250/013	250/015	012/010	012/015	012/016	012/018	226/023	019/010	019/014	019/016	032/013	032/015	031/035	001/010	001/013	001/015	001/018	002/016	237/009
L(mm)	4	7.5	7.5	7.5	10	10	1	1.5	1.4	1.7	7	2.2	3.5	3.5	1.5	1.7	5.7	0.9	1.2	1.3	1.6	3.5	2
Item Code (SHDP+)	462R	364R	265R	165R	144	330	230	130	325	111	246	117	319	145	148	149	248	121					
Short																							
SuperShort																							
Robot Points														0893				0896		0897			
ISO#	237/011	237/015	237/021	237/025	237/032	041/025	041/043	041/050	304/028	068/045	162/033	254/032	243/037	277/029	039/018	038/039	005/018	552/026					
L(mm)	1.7	4.5	5	6.5	5	0.5	0.6	0.7	0.8	1.7	4	5	3	4.9	5	4	7.5	3.6					



Diamond Points FG  
Fine Grit



Designed utilising ne diamond particles to ensure a smooth uniform surface.  
Available in: 12 different shapes with some shapes in short shank (categorised under 1000 class).

Item Code (SHDP+)	F201R	F104R	F102R	F105R	F204	F104	F215	F216	F440R	F414SS	F414S	F414
Short	F1201R				F1204		F1215					
SuperShort												
Robot Points												
ISO#	222/015	193/015	223/016	223/021	165/012	166/015	249/013	249/016	138/007	150/010	150/012	150/014
L(mm)	7	10	9	11.5	8	10	7.5	7.5	0.7	9	9	9

Single \$8.70 10 Pack 40% OFF 100 Pack\* 50% OFF



### Diamond Points FG Super Fine Grit



This range of diamond points have utilised ultra-fine diamond particles to ensure a smooth surface. **Available in: 30 different shapes**, some in short shank (categorised under 1000 class).

Item Code (SHDP+)	SF201	SF101	SF301CR	SF202CR	SF201CR	SF101CR	SF103CR	SF206CR	SF207CR	SF104R	SF102R	SF102RD	SF106RD	SF107RD	SF201RA	SF206RA
Short								825V	826V	SF1104R		833V	836V	839V		
SuperShort																
Robot Points																
ISO#	171/015	173/015	544/018	544/021	545/018	553/019	553/023	544/027	544/031	193/015	223/016	198/018	198/021	198/023	218/016	218/012
L(mm)	7	10	4	5	7	10	10	3.5	3.5	10	9	9	9	9	7	7

Item Code (SHDP+)	SF204	SF215	SF217	SF114	SF318	SF319	SF219	SF416	SF150	SF151	SF440	SF265R	SF313R	SF145
Short				808V					899V					893V
SuperShort														
Robot Points														
ISO#	165/010	249/012	249/013	158/013	254/023	243/014	243/019	249/010	031/035	028	001/015	237/020	284/016	277/029
L(mm)	8	7.5	6	8	2	3.4	5	4	5.7		1.3	5	3.3	4.9

Single \$8.70

10 Pack 40% OFF

100 Pack\* 50% OFF

\* Same Bur



### Diamond Points FG Coarse Grit



Designed utilising larger diamond particles compared to regular grit. **Available in: 8 different shapes.**

Item Code (SHDP+)	C101	C202	C102R	C104LL	C165R	C101X	C103X	C105RX
Short								
SuperShort								
Robot Points								
ISO#	173/017	170/022	222/018	160/022	237/025	173/017	173/023	223/023
L(mm)	10	5	9	10	7	10	10	12

Single \$11.20

10 Pack 40% OFF

\* Same Bur



### Diamond Points FG Super Coarse Grit



These diamond points are designed with large diamond particles. **Available in 10 different shapes.**


Item Code (SHDP+)	SC101	SC102R	SC104R	SC106RD	SC201R	SC202	SC301	SC111	SC144	SC165R
Short										
SuperShort										
Robot Points										
ISO#	018	019	018	024	018	023	018	047	034	026
L(mm)	10	9	10	9	7	5	4.5	1.7	5	6.5

Single \$12.20

10 Pack 40% OFF

\* Same Bur

## SHOFU Flat Fissure Friction Grip Burs



Popular!

Equivalent to the FG541 or ISO code: 109.010 Flat Fissure Style

- Offers optimum material reduction
- Does not clog or generate heat
- Made with high grade diamond particles
- Manufactured with a high level of quality control

100 per bag (bulk pack)  
SD-541

**\$385.60**

## SHOFU Crown & Bridge Preparation Kit



An all you need kit for speedy and accurate tooth preparation for a wide range of crown and bridge restorations, veneers and onlays.

Crown & Bridge Preparation Kit  
SH-F270

**\$158.85**

Usually \$176.50

## SHOFU All Ceramic Preparation Kit



For preparation of all ceramic restorations. 17 Exceptional Robot diamonds.

- Use for anterior, premolar and molar preparations
- Two grits of Robot diamonds: regular and super fine
- Also for inlay and laminate veneer preparations

All Ceramic Preparation Kit  
SH-0955

**\$178.10**

Usually \$197.90

## SHOFU Direct Cosmetic Restoration Kit



A carefully selected assortment of diamond points to help achieve accuracy and efficiency during tooth preparation and contouring / finishing of direct cosmetic restorations.

**The kit contains:**  
8 diamond points FG in 3 grits (regular, fine and superfine)

Direct Cosmetic Restoration Kit  
SH-F270

**\$158.80**

Usually \$176.50

# ethOss® Grow Stronger™ Degranulation Burs

### Key Features & Benefits:

- Designed morphologically to remove soft tissues without harming the bone
- Quick and effective debridement
- Precision manufacture – diamond edges do not grab the bone compared to traditional tungsten burs
- Excellent for preparing a defect for bone grafting with ethoss®

### Advantages:

- Optimal debridement of the inflammatory tissue – cannot be achieved manually
- Multi task drills can be used in implantology, periodontal surgery, endodontic surgery, and root removal
- Maximum safety when working near anatomical structure
- Reusable (more than 100 procedures)

### EthOss Degranulation Refill Kits (x3 burs)

Bur Size	Code
1.0mm	ET-GR-REFILL1.0
2.5mm	ET-GR-REFILL2.5
3.0mm	ET-GR-REFILL3.0
3.5mm	ET-GR-REFILL3.5

**\$400.00**

### EthOss Eco Degranulation Kit

4x Burs & Stand  
1, 2.5, 3, 3.5mm

**\$500.00**

ET-ECGRKIT1

### EthOss Ek Degranulation Burs Kit

**\$673.20**

Includes 4 burs, stand, suction tip & cleaning brush

ET-TEKDG





## SHOFU GIC Polishing Kit

The kit contains all important instruments for finishing and polishing restorations from conventional glass ionomer cements, resin modified glass ionomers and compomers.



**The kit contains:**

- T&F Hybrid Points: 1 each 7406, 7803
- Dura-White: 1 each CN1, FL2, RD1, RD2
- CompoSite: 1 each cup, knife, midpoint
- Ceramisté Ultra II: 1 each cup, knife, midpoint

GIC Polishing Kit

SH-0323

**\$79.80**

Usually \$88.70

## SHOFU Composite Finishing Kit

This kit is suitable for finishing direct resin as well as glass ionomer restorations. Designed for contouring and finishing.



**The kit contains:**

- Dura-Green: 3 each CN1, FL2
- Dura-White: 3 each CN1, FL2

CA shank  
SH-0306

**\$53.20**

Usually \$59.20

FG shank  
SH-0302

**\$51.30**

Usually \$57.10

## SHOFU Porcelain Adjustment Kit

Designed for chairside adjustment and glaze like polishing of porcelain restorations prior to placement.



**The kit contains:**

- Dura-White: 1 each CN1, FL3, RE1
- Ceramisté Standard: 1 each KN7, PC2, WH6
- Ceramisté Ultra: 1 each KN7, PC2, WH6
- Ceramisté Ultra II: 1 each KN7, PC2, WH6

Porcelain Adjustment Kit

SH-0301

**\$51.30**

Usually \$57.10

## SHOFU Composite Polishing Kit

The kit contains all instruments for contouring, finishing and polishing all types of composite restorations and glass ionomer: Dura-Green stones to reduce excess material, Dura-White stones for finishing, CompoSite polishers for prepolishing and CompoSite Fine instruments for high gloss without damage to the enamel.



**The kit contains:**

- Dura-Green: 1 each CN1, FL2
- Dura-White: 1 each CN1, FL2
- CompoSite: 2 each bullet, knife, wheel
- CompoSite Fine: 1 each cup, disk

Composite Polishing Kit

SH-0310

**\$51.30**

Usually \$57.10

## SHOFU Porcelain Veneer Kit

The Porcelain Veneer Kit contains all necessary rotary instruments for the preparation of the tooth, for fine adjustments to porcelain and the final super polish to a high gloss.



**The kit contains:**

- Robot Points, standard: 1x 0835-1
- Robot Points, fine: 1 each 835F-1, 883F-1
- T&F Hybrid Points: 1 each 7404, 7803, 7901
- Dura-White: 1 each FL2, RE1
- Ceramisté Standard: 1 each knife, wheel
- Ceramisté Ultra: 1 each knife, wheel, cup
- Ceramisté Ultra II: 1 each knife, wheel, cup

Porcelain Veneer Kit


SH-0300

**\$51.70**

Usually \$57.50

## SHOFU CAD/CAM Bur Preparation Kit\*

Preparation kit for all-ceramic and hybrid ceramic CAD/CAM restorations.



Eleven diamond coated and shape adjusted rotary instruments for the quick and easy conservative preparation as well as for the preparation for all-ceramic or hybrid ceramic restorations; anterior and posterior. The high-quality bur stand is made of stainless steel and can be sterilised.

CAD/CAM Kit  
SH-079CAD

**\$47.00**

Usually \$52.00

\*on backorder



### Bur Station II

Autoclavable and economical plastic bur blocks.



- High thermal resistance
- High durability against chemicals
- Holds 11 FG shank and 6 CA shank burs

SH-0797 Blue  
SH-0796 Yellow

**\$22.90**  
Usually \$25.50

### Bur Station

A handy and convenient bur stand specially designed for easy sterilisation and storage of FG burs.

Sterilised by autoclaving to 135 °C (270 °F), ultrasonic cleaning and cold sterilisation.

Unique design can hold 12 FG burs  
Available in: Blue

Bur Station  
SH-0798

**\$17.40**  
Usually \$19.40

## DIRECTA Bur Stand Roto




- Easy replacement of burs or files
- Safe handling - minimal risk of injury
- Rotates to facilitate quick and easy selection
- Close fitting lid to protect burs from dust and contamination
- Cold sterilisation

DA-4122-85 FG - Holds up to 60 burs  
DA-4122-87 Holds 18 FG & 44 HP/RA burs


**\$50.40**  
Usually \$56.00




### Aluminium Bur Stands



24 FG BURS 60x40mm  
ME-965-1




24 RA BURS 60x40mm  
ME-965-2




12 FG / 12 RA BURS 60x40mm  
ME-965-3


**\$32.10**  
Usually \$35.70



36 FG BURS 60x40mm  
ME-966-1




36 RA BURS 60x40mm  
ME-966-2




18 FG / 18 RA BURS 60x40mm  
ME-966-3


**\$36.70**  
Usually \$40.80



48 RA Burs 75 x 40mm  
ME-967-1



48 RA Burs 75 x 40mm  
ME-967-2



24FG / 24RA 75 x 40mm  
ME-967-3

**\$22.90**  
Usually \$25.50

# DIRECTA

## Bur Box

- 48 x 27 x 9mm
- Aluminium
- Autoclavable



**\$24.30**

DA-4160-84

Usually \$27.00

## Bur Stand

- 78 x 40 x 7mm
- Holds 24 CA and FG burs
- Autoclavable



**\$45.00**

DA-4163-05

Usually \$50.00

## Bur Stand Mini

- 12 HP / RA-12 FG
- Autoclavable



**\$62.10**

DA-4163-07

Usually \$69.00

# DIRECTA

## PractiPal®

Small file and bur holders in two different models to minimise sharps injuries. All bur and file holders may be used with PractiPal Full and Half Tray or as a stand alone and can easily be secured onto the PractiPal Tray through the whole cleaning and disinfection process. The file stand holds 13 files while the bur stand holds 2 long burs, 5 RA and 8 FG burs.



## Compact Bur Stand Mint



**\$45.00**

DA-1151-23

Usually \$50.00

## Compact File Holder Blue



**\$45.00**

DA-1151-83

Usually \$50.00

**Additional colours available on request**

See page 30 for PractiPal® Trays...

INTRODUCTORY OFFER!



Saves Time



Easy to Use



Versatile



Reduces Risk of Cross Infection



Patient Comfort



## Master the art of precision with Pro-Matrix

Pro-Matrix is a single-use solution for amalgam and composite restorations for use in all quadrants in every class. Available as a standard matrix band and a contoured band, one user's comment sum Pro-Matrix up...

"These matrix bands are pre-assembled. They save time and are very convenient."

- Fully disposable
- Single Product System no assembly or application tools required
- Circumferential Band can be used for MO/OD/MOD/ Cusp replacements
- High-tension mechanism for better compacting
- Easy contouring to tooth shape
- Narrow front section better visibility and patient comfort
- Suitable for all quadrants
- Suitable for use with all restorative materials
- Suitable for use with wedges and rubber dams



"Pro-Matrix is my preferred matrix system due to ease of use, versatility and excellent clinical outcome."

**Ash Rayarel BDS MFGDP(UK) Msc.**

Ash is the Principal Dentist and owner of Ash Dental Practice and Ash Aesthetics.



AK-19000 NARROW

AK-19005 WIDE

Normally \$127.50 **\$114.70** ea




Pro-Matrix Curve has the additional benefits of:

- Pre-contoured to tooth shape
- Longer band for larger teeth

AK-19001 CURVE NARROW

AK-19006 CURVE WIDE

Normally \$146.90 **\$132.20** ea



**Clinician's Choice**  
Designed for better dentistry.

**Convexi-T™**  
Convex Tofflemire® Matrix Bands

**Convexi-T S2™**  
Stainless Steel Convex Tofflemire® Matrix Bands


### Produce Naturally Contoured Interproximal Contacts

Using regular matrix bands to restore interproximal restorations is always difficult because most matrix bands are flat occluso-gingivally. Even when the contact is burnished, the matrix band seldom holds its shape. This leads to flat interproximal contours, where the contact point is translocated high occlusally. This contact is often broken when creating occlusal embrasure form, or prematurely breaks down under normal occlusal forces.

**Convexi-T™ Convex Tofflemire™ Matrix Bands** are pre-contoured buccal/lingually and occlusal/gingivally and are available in two options. **Convexi-T** matrix bands are 0.025mm thick and are dead soft, allowing them to be further contoured if necessary. **The Convexi-T S2 (stainless steel)** matrix bands are slightly more rigid at 0.04mm thick, facilitating easy placement when preparations are tight - even having the ability to push through a partially closed preparation. Both options will provide a tighter and more naturally contoured contact for your Class II restorations.


Both versions of **Convexi-T** matrix bands, available in widths of 5.5mm and 6.3mm, facilitate the creation of interproximal contours necessary to recreate the contact areas and natural convex anatomy of the interproximal area.

**\$226.80**  
Usually \$252.00



Convexi-T Matrix Band  
6.3mm X 0.025mm (4X30pk)  
CC-206105

Convexi-T Matrix Band  
5.5mm X 0.025mm (4X30pk)  
CC-206106



Convexi-T S2 Matrix Band  
6.3mm X 0.04mm (4X30pk)  
CC-206107

Convexi-T S2 Matrix Band  
5.5mm X 0.04mm (4X30pk)  
CC-206108



## Matrix Band Holders

**Tofflemire Matrix Band Holder**

**Universal**  
5-6-7mm

ME-4975

**\$45.90**  
Usually \$51.00




**Tofflemire Matrix Band Holder**

**Junior**  
5mm

ME-4976

**\$45.90**  
Usually \$51.00




**Siqveland Matrix Retainer**

**Wide**  
5-6-7mm

ME-4972-WL

**\$50.40**  
Usually \$56.10



**Siqveland Matrix Retainer**

**Narrow**  
5-6-7mm

ME-4972-NS

**\$50.40**  
Usually \$56.10



## Tofflemire Matrix Bands

**0.05mm 12 Pack**

ME-4979-12

**\$13.70**  
Usually \$15.30



**0.05mm 12 Pack**

ME-4979-13

**\$13.70**  
Usually \$15.30



## BIOCLEAR Diamond Wedge

The Diamond Wedge requires less effort to place than traditional wedges due to the diamond shaped cut-out in the tip, which collapses during insertion. The wedge then springs open once through the embrasure, creating a tight seal and preventing the wedge from backing out.

- Easy to place
- Will not back out
- Strong separation pressure
- Atraumatic to papilla




**\$87.30**  
Usually \$97.00

50 Packs

BC-K-DW  
**Assorted Kit 80pcs**  
16 of each size

**\$124.20**  
Usually \$138.00

Small Pink	BC-DW-S-50
Medium Orange	BC-DW-M-50
Large Yellow	BC-DW-L-50
Extra Large Blue	BC-DW-XL-50
Deep Caries Green	BC-DW-DC-50



## DIRECTA Barman's Wedges

Barman's Anatomical Wedges achieve better and safer procedures thanks to

its anatomical shape. The role of the wedges is to attach the matrices to their correct position and to achieve a better separation from the adjacent tooth.

The wedge is progressively thicker from pointed end making the wedge fit most spaces. Available in two sizes small (natural) and large (green). Made of birch wood which makes the wedge soft and pliable.



Barman's Easy Grip	DA-6027-90	
Barman's Anatomical Mixed	DA-6027-91	<b>200 Packs</b>
Barman's Anatomical Small	DA-6027-92	<b>\$37.80</b>
Barman's Anatomical Large	DA-6027-93	Usually \$42.00

## Contoured Wooden Wedges



Clinician's Choice® Contoured Wood Wedges have precise anatomical contour and can help you develop great posterior contacts.

- Made from maple, for strength and compressibility
- Prevents dislodgement or movement of the matrix band during composite placement and packing
- Available in 11mm and 12mm, the smallest wood wedge sizes available in dentistry. Also available in 13mm, 15mm, and 17mm
- Can be trimmed without splintering for a custom fit
- Wedges' colour coding agent will not discolour the tooth or restoration



CC-095940  
**Assorted Pack 400pcs**

150 x 11mm (Pink)  
100 x 12mm (Blue)  
100 x 13mm (Green)  
50 x 15mm (Yellow)

**\$135.90**  
Usually \$151.00

400 Pack Refills

11mm Pink	CC-095411
12mm Blue	CC-095412
13mm Green	CC-095413
15mm Yellow	CC-095415
17mm White	CC-095417

**\$78.30**  
Usually \$87.00

## DualForce Active-Wedges



DualForce™ Active-Wedges actively engage the matrix band and adjacent tooth.

Active-Wedges are designed to easily glide into position with a collapsing tip that re-opens once fully placed.

The unique Active-Wedge contour and strong internal spine places significant force on the matrix band to seal margins, eliminating overhangs and preventing composite flash - translating to less finishing time.



CC-209000  
**Assorted Pack 200pcs**  
40 of each size

**\$260.10**  
Usually \$289.00

100 Packs

Size 1 Purple	CC-209001
Size 2 White	CC-209002
Size 3 Blue	CC-209003
Size 4 Green	CC-209004
Deep Seal Black	CC-209005

**\$171.00**  
Usually \$190.00

## DIRECTA FenderWedge®

Protects the adjacent tooth & gingiva during preparations.

- Confident preparations
- Protects the adjacent tooth
- Pre-Separation
- Protects gingiva




Bur damage to adjacent teeth is a common problem in everyday dental practice. Research shows that 70% of adjacent teeth suffer damage during class II preparations.\* FenderWedge prevents damage to the adjacent tooth and gingiva that commonly

occurs during Class II preparations.


**FenderWedge** is a combination of a steel plate and a plastic wedge. It facilitates easy application of a matrix. To avoid aspiration it is recommended to secure FenderWedge with a waxed dental floss.




**Extra Small**  
36 Pack DA-6027-54  
100 Pack DA-6028-03



**Small**  
36 Pack DA-6027-51  
100 Pack DA-6028-00



**Medium**  
36 Pack DA-6027-52  
100 Pack DA-6028-01



**Large**  
36 Pack DA-6027-53  
100 Pack DA-6028-02

DA-6027-50  
**Assorted Pack 144pcs**  
36 of each size  
**\$194.40**  
Usually \$243.00



**36 Pack \$56.80**  
Usually \$71.00

**100 Pack \$119.20**  
Usually \$149.00

\* Qvist V, Johannessen L, Bruun M (J Dent Res 1992 Jul; 71 (7): 1370-3) Lussi A, Gygax M. (J Dent. 1998 Jul-Aug; 26 (5-6):435-41), Medeiros VA, Seddon RP (J Dent. 2000 Feb; 28 (2): 103-10)

## DIRECTA FenderMate®


Sectional matrix for Class II composite restorations.

- Predictable contact and anatomic shape
- One piece design
- No cervical overhang
- Designed for minimally invasive dentistry




**FenderMate** is a pre-curved, one-piece sectional matrix and wedge that provides for quick, safe and predictable composite restorations, with a tight contact and cervical sealing. FenderMate is designed to be inserted either buccally or lingually.


The matrix reaches from the base of the wedge to just above the occlusal surface. The side of the wedge facing the adjacent tooth has an angled wing.



**Left Regular**  
18 Pack DA-6027-71  
100 Pack DA-6028-04



**Left Narrow**  
18 Pack DA-6027-72  
100 Pack DA-6028-05



**Right Narrow**  
18 Pack DA-6027-74  
100 Pack DA-6028-07



**Right Regular**  
18 Pack DA-6027-73  
100 Pack DA-6028-06

DA-6027-70  
**Assorted Pack 72pcs**  
18 of each size  
**\$304.00**  
Usually \$380.00



**18 Pack \$87.20**  
Usually \$109.00

**100 Pack \$374.40**  
Usually \$468.00


## DIRECTA FenderMate® Prime

Protective and easy matrix for primary teeth.

- Designed for primary teeth
- Easy placement
- Bendable
- One piece design

**FenderMate Prime** uniquely combines a preparation shield and matrix for use in Class II fillings of primary teeth, where quick and simple solutions are highly valued. It is composed of an interproximal plastic wedge with an inbuilt steel plate. The "boat-shaped" tip allows the gingiva to be

compressed without causing damage and reduces the risk of bleeding. FenderMate Prime is available in two sizes; long and short. FenderMate Prime green (long) is bendable and can be shaped to the contour of the tooth. FenderMate Prime yellow has a shorter, more rigid plate for very tight interproximal spaces.



**Prime Long**  
18 Pack DA-6028-17  
100 Pack DA-6028-09



**Prime Short**  
18 Pack DA-6028-16  
100 Pack DA-6028-08

DA-6028-15  
**Assorted Pack 36pcs**  
18 of each size  
**\$87.20**  
Usually \$109.00



**18 Pack \$40.00**  
Usually \$50.00

**100 Pack \$132.80**  
Usually \$166.00

# BT MATRIX



## Black Triangle Matrix System

Black Triangle matrices are custom designed for black triangle treatment and require no trimming prior to seating. A colour-coded gauge and matrices make choosing the correct curvature for gingival embrasure closure simple and predictable. Curvature is indicated by colours at the top of the matrix and correspond with a gauge that, when inserted into the black triangle, indicates which curvature to use.



**Intro Kit**  
**\$673.00**

**BC-801003** Usually \$748.00



Item	Quantity
Black Triangle Matrix – Small Blue	10
Black Triangle Matrix – Small Green	10
Black Triangle Matrix – Small Yellow	10
Black Triangle Matrix – Small Pink	10
Black Triangle Matrix – Large Blue	10
Black Triangle Matrix – Large Green	10
Black Triangle Matrix – Large Yellow	10
Black Triangle Matrix – Large Pink	10
Black Triangle Gauge	1
TruContact Heavy Sander, Double-Sided, Red	1
TruContact Light Sander, Single-Sided, Orange	1
TruContact Saw, Blue	1
Bioclear Dual Colour Disclosing Solution – 3ml	1

**\$198.00**

**Procedure Kit** **BC-801010** Usually \$220.00

Item	Quantity
Black Triangle Matrix - Small Blue	3
Black Triangle Matrix - Small Green	3
Black Triangle Matrix - Small Yellow	3
Black Triangle Matrix - Small Pink	3
Black Triangle Matrix - Large Blue	3
Black Triangle Matrix - Large Green	3
Black Triangle Matrix - Large Yellow	3
Black Triangle Matrix - Large Pink	3
Black Triangle Gauge	1

### Complete HD Anterior Matrix Kit

Item	Quantity
HD Matrix A Series – A101	15
HD Matrix A Series – A102	15
HD Matrix A Series – A103	15
HD Matrix A Series – A104	15
HD Matrix A Series – A105	15
HD Matrix DC Series – DC201	15
HD Matrix DC Series – DC202	15
HD Matrix DC Series – DC203	15
Diamond Wedge – Small	8
Diamond Wedge – Medium	8



**BC-820197**

**\$613.80**

Usually \$682.00

### HD Anterior Matrix Starter Kit

Item	Quantity
HD Matrix A Series – A101	15
HD Matrix A Series – A102	15
HD Matrix A Series – A103	15
HD Matrix DC Series – DC201	15
HD Matrix DC Series – DC202	15
HD Matrix DC Series – DC203	15
Diamond Wedge – Small	8
Diamond Wedge – Medium	8



**BC-820196**

**\$544.00**

Usually \$605.00



# BIOCLEAR

## Posterior Matrix Kits

### Biofit HD Posterior Complete Kit

Kit contains the entire range of Biofit HD Matrices (appropriate for 75% of posterior molar cases) and also includes Diamond Wedges (which have strong separation yet require less effort to place than traditional wedges) and Twin Rings (comprised of strong, long-lived twin NiTi wires).



BC-818081

**\$881.00**

 Usually \$979.00 **Kit Contains:**

- |  |  |
|--|--|
| 25x Posterior Matrix Series - Molar 4.5    | 16x Diamond Wedge Series - Extra Large |
| 25x Posterior Matrix Series - Molar 5.5    | 16x Diamond Wedge Series - Deep Caries |
| 25x Posterior Matrix Series - Molar 6.5    | 1x Twin Ring Series - Molar            |
| 25x Posterior Matrix Series - Premolar 5.5 | 1x Twin Ring Series - Premolar         |
| 16x Diamond Wedge Series - Small           | 1x Twin Ring Forceps                   |
| 16x Diamond Wedge Series - Medium          | 1x Curved Micro Pliers                 |
| 16x Diamond Wedge Series - Large           |  |

### Biofit HD Posterior Mini Assorted Kit


**Kit Contains:**

- 25x Posterior Matrix Series - Molar 4.5
- 25x Posterior Matrix Series - Molar 5.5
- 25x Posterior Matrix Series - Molar 6.5

BC-818001

**\$247.50**

Usually \$275.00

### Biofit EVOLVE Posterior Kit


**Kit Contains:**

- 10x Bioclear Evolve Blue 5mm
- 10x Bioclear Evolve Blue 6mm
- 10x Bioclear Evolve Blue 7mm
- 10x Bioclear Evolve Blue 8mm
- 10x Bioclear Evolve Blue 9mm
- 10x Bioclear Evolve Orange 6mm
- 10x Bioclear Evolve Orange 7mm
- 10x Bioclear Evolve Orange 9mm
- 10x Bioclear Evolve Pink 6mm
- 10x Bioclear Evolve Pink 7mm
- 10x Bioclear Evolve Pink 9mm

BC-827061

**\$379.80**

Usually \$422.00

### Biofit EVOLVE Posterior Complete Kit

The Kit offers everything you need to start using Bioclear Evolve Matrix Series in your practice. Following the legacy of the Biofit HD Matrix Series, the Evolve Matrices have the same benefits and additional features for posterior cases.


**\$1,140.30**

BC-827051

**Kit Contains:**

- |                                   |  |
|-----------------------------------|--|
| 25x Bioclear Evolve Blue 5mm      | 16x Diamond Wedge Series - Large           |
| 25x Bioclear Evolve Blue 6mm      | 16x Diamond Wedge Series - Extra Large     |
| 25x Bioclear Evolve Blue 7mm      | 16x Diamond Wedge Series - Deep Caries     |
| 25x Bioclear Evolve Blue 9mm      | 1x Twin Ring Series - Molar                |
| 25x Bioclear Evolve Orange 7mm    | 1x Twin Ring Series - Premolar             |
| 25x Bioclear Evolve Orange 9mm    | 1x Twin Ring Forceps                       |
| 16x Diamond Wedge Series - Small  | 1x Bioclear Push-Pull                      |
| 16x Diamond Wedge Series - Medium | 1x Bioclear Evolve Matrix Height Indicator |

**MEDESYS** ITALY Rubber Dam Punches

Rubber Dam Forceps Ivory




ME-5551  
**\$82.60**  
 Usually \$91.80

Rubber Dam Punch Ainsworth



ME-5560-2  
**\$160.65**  
 Usually \$178.50

**MEDESYS** ITALY Rubber Dam Frames



- ME-5554 Small 80 x 90mm (Young)
- ME-5553 Medium 110 x 110mm (Young)
- ME-5552 Large 140 x 140mm (Young)

All Sizes  
**\$36.70**  
 Usually \$40.80

**MEDESYS** ITALY

Rubber Dam Medium – Mint



High quality low powder latex dental dams. Processed to low protein content of <50 micrograms/dm<sup>2</sup> of material.

- ME-5610-M-V Green Medium (36 Sheets)
- ME-5610-M-N Black Medium (36 Sheets)

**\$33.00**  
 Usually \$36.70

**MEDESYS** ITALY Dental Dam Kit



GAMMAFIX trays by MEDESYS feature high quality and functionality; they are specifically designed to facilitate the handling of instruments with maximum safety and to guarantee sterilisation that is safe for the patient but does not damage the instrument itself.

Dental Dam Kit Includes:

981/20	Gammatrix	5595/12A	Clamp No. 12A
	Cinc tray	5595/13A	Clamp No. 13A
5550	Dam punch	5595/27	Clamp No. 27
5551	Clamp forceps	5595/29	Clamp No. 29
5552	Dam frame	5595/210	Clamp No. 210
5553	Dam frame	5595/212	Clamp No. 212
5554	Dam frame	5595/W8A	Clamp No. W8A
5595/00	Clamp No. 00	5605	Clamp tray
5595/3	Clamp No. 3	5610/M-V	Rubber dam box (not shown)
5595/4	Clamp No. 4		

**\$872.10**  
 Usually \$969.00

ME-1675-10

**centrix**  
MAKING DENTISTRY EASIER

# Exposé



## Caries Indicators

Exposé is the simplest-to-use caries detector product. Instead of messy squeeze bottles or syringes, Exposé is pre-dried on a Benda® Micro applicator. Simply touch it to a wet tooth surface to apply. No drip or mess that can stain clothing or soft tissue.

- Perfect for minimally invasive preps; allows you to identify and remove affected tooth structure, leaving healthy enamel and dentine
- Single-use applicators make tray set-up easy – simply place one on the tray – no syringe, no sheath, no tips to replace, no cross-contamination
- Contains an FD&C dye in a propylene glycol base
- Works by staining the by-products of carious lesions

**\$110.10**  
Usually \$122.40

**CE-310334** Blue Fine  
**CE-310335** Blue Extra Fine

MAKING DENTISTRY EASIER®

**VISTA** | APEX

# Vista-Green™ Caries Indicator

Vista-Green™ Glycol-based solution to facilitate removal of carious dentine. Ensures maximum amount of inner dentine is preserved and infected dentine is removed.

Easily locate root canal orifices and calcified canals.

**\$85.50**  
Usually \$95.00

**V-502750**

**Standard Kit includes:**

- 4x 1.2mL prefilled syringes
- 20x Sol-u-Flo™ tips



**V-502755**

**Refill Kit includes:**

- 4x 1.2mL prefilled syringes

**\$76.50**  
Usually \$85.00

# PULPDENT®

## SNOOP™ Caries Detecting Dye

For Conservative Cavity Preps

Snoop™ distinguishes between outer infected and inner affected dentine in 10 seconds.

**REALITY**  
Four Star Award



Snoop stains infected carious dentine in 10 seconds and helps preserve vital dentine that should not be removed. This is an important tool for conservative dentistry.

- Detects caries in difficult to see places
- Dark blue colour provides strong contrast with dentine and the pulp
- Propylene glycol formula is reported to be more accurate than water-based formulas
- Economical dropper bottle
- Rated Four Stars by REALITY

**\$76.50**  
Usually \$85.00

12ml Bottle  
**PL-DSNOOP**

# PULPDENT®

## DENTINE DESENSITISER

Contains 5% Glutaraldehyde, Fluoride and Water. Apply to all dentine surfaces, including cavity, crown and inlay preparations and cervical areas. Compatible with all adhesives, composites, provisional and permanent dental cements.



Pulpdent Dentine Desensitiser contains 5% glutaraldehyde in water. It is not a dentine primer; however, it can be used in conjunction with adhesives and composites, as well as with provisional and permanent traditional dental cements. It does not interfere with bonding techniques and does not affect bond strengths of adhesives, composites and resin cements. Fluoride is added to the Desensitiser to enhance the stability of the solution.

12ml Bottle  
**PL-DEES**

**\$91.80**  
Usually \$102.00

## Clinician's Choice

### Max Etch 35% Phosphoric Acid

Max Etch has a 35% concentration of phosphoric acid, which represents the gold standard in achieving a consistent, micro-mechanically retentive enamel and dentine surface for virtually all bonding procedures.



Max Etch  
35% 5ml  
CC-345503  
**\$26.00**  
Usually \$29.00

Max Etch  
35% Kit  
30ml  
CC-345302  
**\$62.00**  
Usually \$69.00

## PULPDENT®

### Etch Rite 38% Phosphoric Acid

Etch-Rite is a superior quality dental etching gel specially formulated for maximum working and handling characteristics.

- Thixotropic
- Stays where placed - will not run onto exposed tissue
- Washes off easily without leaving any residue
- 20 second etch for enamel and 15 second etch for dentine



Etch Rite  
4x1.2ml  
PL-D-ETCH  
**\$28.80**  
Usually \$32.00

Etch Rite  
Jumbo 2x  
25ml  
PL-DET50R  
**\$63.00**  
Usually \$70.00

Etch Rite  
Bulk Pack  
24x 1.2ml  
PL-DET-24  
**\$128.70**  
Usually \$143.00



### Best-Etch 37% Phosphoric Acid

Best-Etch™ is a 37% Phosphoric Acid Etch, designed to etch and prepare dentine and enamel.

- Viscous solution that stays in place
- Rinses clean with water and leaves no residue
- Pre-bent applicator tips ensure precise delivery



Best Etch  
37 % 30ml  
V-502230  
**\$92.70**  
Usually \$103.00

Best Etch  
37 % Bulk  
30ml  
V-502235  
**\$106.20**  
Usually \$118.00

Best Etch  
37 %  
Value Pack  
12 x 3ml  
V-502200  
**\$106.20**  
Usually \$118.00

## PULPDENT®

### Etch Royale 37% Phosphoric Acid

For clinicians who prefer a creamier gel that readily settles into dentine and enamel, but does not run, Etch Royale is the perfect choice. Etch Royale has all the same features as Etch-Rite™, but the consistency is slightly creamier than its famous sister product.



Etch Royale  
4x1.2ml  
PL-DER  
**\$36.00**  
Usually \$40.00

Etch Royale  
Jumbo Refill  
2x25ml  
PL-DER50R  
**\$63.00**  
Usually \$70.00



### Porcelain Etch 9.5% Hydrofluoric Acid

Viscous 9.5% Hydrofluoric Acid Solution. For porcelain veneers and inlays prior to cementation or repairs.

- Viscous solution
- Rinses clean away
- Includes adjustable Spira-Flo™ brush tips with curved cannula for optimal entry
- Silane solution used after etchant enhances the bonding of resin material



Porcelain  
Etch 9.5%  
2x 1.2ml  
V-502500  
**\$81.00**  
Usually \$90.00

## PULPDENT®

### Porcelain Etch Gel 9.6% Hydrofluoric Acid

A one-minute application etches porcelain crowns, veneers, inlays and fractured porcelain surfaces prior to bonding. This superior quality porcelain etch gel will not stain porcelain, composites or resin cements.



Porcelain  
Etch 4x1.2ml  
PL-DPEG  
**\$93.60**  
Usually \$104.00



MPa Max is a total-etch adhesive system that provides maximum bond strength and maximum prevention of post-operative sensitivity. MPa Max formulation includes 0.2% chlorhexidine to improve bond longevity; inhibiting MMPs. MPa Max's simple technique requires a single application on each adhesive wall.

5ml Refill Bottle

**\$161.00**

Usually ~~\$179.00~~

CC-266506

# MPa™ Max

Maximum Performance Adhesive

1 x 5mL bottle MPa Max, 1 x 5mL syringe Max Etch, 1 x 5mL bottle G5 All-Purpose Desensitiser, tips, mixing wells, applicator brushes, instructions/SDS

MPa Max Starter Kit: **\$244.80**

CC-266901 Usually ~~\$272.00~~



**Clinician's Choice**  
Designed for better dentistry.

## BeutiBond Adhesive Kit

Offers a simple, safe and efficient application technique with high initial bond strength comparable with two step bonding systems

HEMA-free composition and unique dual adhesive monomers ensures secure bond to enamel and dentine while minimising post-operative sensitivity with no gingival blanching and preventing aging caused by hydrolysis of the hybrid and adhesive layers

A very thin, bubble-free, 5 µm bonding layer helps achieve superior aesthetics with a more conservative approach in flat and shallow cavities with a secure marginal seal

Short and convenient application time of 30 seconds



**30% OFF**

6 ml Bottle, V-Dish x 25, Microtip Fine Handle x 1, Application Tips x 50

**\$99.20**

Usually ~~\$141.80~~

SH-1781



## Latest Product! BeutiBond Xtreme

The True Universal Adhesive  
Etch, prime and bond in one step!



- For a wide range of applications and adherends with one bottle
- Bonds equally well to wet and dry dentine
- No filler and no HEMA - protection against hydrolysis
- Compatible with self-etch, total-etch and selective-etch techniques
- Simplified procedure for any adherend surface
- No additional primer needed

5ml Bottle, Microbrush x50, V-Dish x25

**\$128.50**

Usually ~~\$142.80~~

SH-Y2448

## FL-BOND II Adhesive Kit

The 2-step, self etching, fluoride releasing adhesive system features a unique Primer and Bonding Agent to provide an excellent bond to both enamel and dentine with a secure marginal seal.

The FL-BOND II PRIMER contains a new effective adhesion promoting monomer and is acetone free with no incorporation of HEMA to minimise odour and post operative sensitivity.

The FL-BOND II BONDING AGENT contains the S-PRG filler (Surface Pre-Reacted Glass-ionomer), which helps to reinforce the bonding interface between the restorative material and the tooth structure.



Bond 5ml, Primer 5ml, Red Microbrushes x25, Yellow Microbrushes x25, Mixing Dish

**\$130.30**

Usually ~~\$144.80~~

SH-1306



# Have you met the

Universal		Enamel / Gingiva		Injectable	
Beautiful II	Beautiful II LS	Beautiful II - Enamel	Beautiful II - Gingiva	Beautiful Injectable X	Beautiful Injectable XSL
					
<p>Universal fluoride releasing aesthetic restorative material.</p> <p>Low Shrinkage Universal fluoride releasing aesthetic restorative material.</p>	<p>4 shades with varying levels of brightness to mimic natural enamel.</p> <p>smooth, creamy consistency and long working time to easily recreate the surface textures of gingival tissue.</p>	<p>Beautiful Injectable X is ideal for restoration of all classes of cavities, minimally invasive cosmetic restorations,</p> <p>Beautiful Injectable X SL - Self-Lavelling</p>			
<p>A1 SH-1401</p> <p>A2 SH-1402</p> <p>A3 SH-1403</p> <p>A3.5 SH-1404</p> <p>A4 SH-1405</p> <p>B1 SH-1408</p> <p>B2 SH-1409</p> <p>B3 SH-1410</p> <p>C2 SH-1412</p> <p>C3 SH-1413</p> <p>A20 SH-1419</p> <p>A30 SH-1420</p>	<p>A1 SH-Y2270</p> <p>A2 SH-Y2271</p> <p>A3 SH-Y2272</p> <p>A3.5 SH-Y2273</p> <p>A4 SH-Y2274</p> <p>B1 SH-Y2275</p> <p>B3 SH-Y2277</p> <p>A20 SH-Y2280</p> <p>A30 SH-Y2281</p> <p>BW SH-Y2283</p>	<p>T SH-Y2250</p> <p>HVT SH-Y2251</p> <p>LVT SH-Y2252</p> <p>AM SH-Y2253</p>	<p>G-LP SH-Y2256</p> <p>G-DP SH-Y2257</p> <p>G-V SH-Y2258</p> <p>G-Br SH-Y2259</p> <p>G-Or SH-Y2260</p>	<p>A1 SH-Y2501</p> <p>A2 SH-Y2502</p> <p>A3 SH-Y2503</p> <p>A10 SH-Y2511</p> <p>A20 SH-Y2512</p> <p>INC SH-Y2513</p> <p>BW SH-Y2514</p>	<p>A1 SH-Y2521</p> <p>A2 SH-Y2522</p> <p>A3 SH-Y2523</p> <p>A3.5 SH-Y2524</p> <p>A30 SH-Y2531</p> <p>INC SH-Y2532</p> <p>BW SH-Y2533</p>
<p>4.5gm Syringe</p> <p>20 0.25 gm Tips</p>	<p>4.0gm Syringe</p> <p>20 0.25 gm Tips</p>	<p>2.5gm Syringe</p>	<p>2.5gm Syringe</p>	<p>2.2gm Syringe</p> <p>5 Needle Tips</p>	<p>2.2gm Syringe</p> <p>5 Needle Tips</p>
<b>\$86.70</b>	<b>\$86.70</b>	<b>\$47.90</b>	<b>\$47.90</b>	<b>\$68.30</b>	<b>\$68.30</b>

# Composite Family?



Flowable		Bulk Fill		Opaquer
Beautiful Flow F02	Beautiful Flow F10	Beautiful Bulk	Beautiful Bulk Flow	Beautiful Opaquer
				
Created to retain the shape and have moderate flow, especially for restoring cervical areas.	Smooth flow, especially for restoring gingival box of class II, Pit & Fissure Sealant and as a cavity liner.	One bulk-fill system – two viscosities for all indications of conservative treatment.	Designed to provide a self-levelling quality for seamless cavity adaptation, and for dentin replacement.	Beautiful Opaquer is used for masking discolourations, staining and metal colour as well as pulp capping.
A1 SH-1431 A2 SH-1432 A3 SH-1433 A3.5 SH-1434 A4 SH-1435 A30 SH-1449 GUM SH-1451	A1 SH-1461 A2 SH-1462 A3 SH-1463 A3.5 SH-1464 A4 SH-1465 A30 SH-1479 A3T SH-1482	Univ SH-2034 A SH-2035	Univ SH-2030 Dent SH-2031	UO SH-1363 LO SH-1364
2.0gm Syringe 5 Needle Tips	2.0gm Syringe 5 Needle Tips	4.5gm Syringe	2.4gm Syringe	2.0gm Syringe 5 Needle Tips
\$63.20	\$63.20	\$88.70	\$88.70	\$75.50


Buy 5  
get 5  
FREE

OR

Buy 2  
get 1  
FREE

# SDI SPECIAL

DECEMBER 2023 - JANUARY 2024



**ionglass™**  
BIOMIMETIC FORMULA



THE RIVA RANGE

**BUY 2 GET 1  
FREE\***

No lock in contracts

Guaranteed savings  
with reduced inventory



**CALL FREE 0508 802 557**

Fax: 04 3873874

[sales@healthcareessentials.nz](mailto:sales@healthcareessentials.nz)

[www.healthcareessentials.nz](http://www.healthcareessentials.nz)

**YOUR LOCAL REPRESENTATIVE:**

Auckland / North Shore / Northland  
Emma-Jane Somerfield 021 818 919

Auckland / Waitakere / Pukekohe  
Fiona Zuccani-Heard 021 922 014

Waikato / Bay of Plenty / Taupo  
Vicki Slow 021 899 729

Wellington / Lower North Island / Gisborne  
Anne Carter 021 450 655

South Island  
Nicky Grant 021 796 360

\*All prices exclude GST. Offers cannot be combined. Valid on invoices dated December 1st - January 31st 2024. All claims for bonus goods must be made within 30 days from the end of promotion date. Bonus goods must be of equal or lesser value. Offers are valid in New Zealand only.



## THE RIVA RANGE

# BUY 2 GET 1 FREE\*

## RIVA SELF CURE/HV

BULK FILL GLASS HYBRID RESTORATIVE SYSTEM

- Contains SDI's proprietary ionglass™, a high ion releasing, radiopaque, bioactive glass that remineralises and prevents secondary caries formation
- Very high fluoride releasing and recharging ability
- Exceptional Handling; Non-sticky, non-slumping and sculptable formulation
- Easy Activation Capsule - 70% easier to activate
- Ideal materials for dentine replacement
- Radiopaque
- BPA & HEMA free

**BOX OF 50 CAPSULES - MRP \$234.60+ GST**



## RIVA LIGHT CURE/HV

LIGHT CURED RESIN REINFORCED GLASS IONOMER RESTORATIVE MATERIAL

- Contains SDI's proprietary ionglass™, a high ion releasing, radiopaque, bioactive glass that remineralises and fights against secondary caries formation
- Exceptional handling; easy to sculpt and non-sticky - available in both regular and high viscosity
- Easy Activation Capsule - 70% easier to activate
- Long term wear resistance due to high surface hardness
- High flexural strengths
- BPA free

**BOX OF 50 CAPSULES - MRP \$292.70 + GST**



## RIVA CEM AUTOMIX

RESIN MODIFIED GLASS IONOMER LUTING CEMENT

- Superior bond (enamel & dentine)
- Highest flexural strength
- More tips / less tip waste
- The best long term aesthetics
- Light tack cure option
- Sustained fluoride release
- BPA free
- Shade: light yellow

**2 X 4G SYRINGE REFILL - MRP \$109 + GST**



## RIVA LUTING

- Chemically adheres to tooth structure and metal substrates
- Very low film thickness - flows easily between restoration and tooth
- Exceptional bond strengths
- BPA & HEMA free

**BOX OF 50 CAPSULES - MRP \$237.70 + GST**



## RIVA LUTING PLUS

- Chemically adheres to tooth and substrates
- Exceptional bond strengths
- High flexural strength - reduced cracking and failure in flexion

**BOX OF 50 CAPSULES - MRP \$268.30 + GST**



## RIVA PROTECT

- ACP nanotechnology - supercharges remineralisation
- Very high fluoride releasing and recharging ability to protect teeth from acid attacks

**BOX OF 50 CAPSULES - MRP \$237.70 + GST**



## RIVA SILVER

MERCURY FREE, FLUORIDE RELEASING, SILVER RESTORATIVE

- Silver reinforced
- Ideal for abfraction lesion restorations
- Ideal for bulk filling
- No shrinkage

**BOX OF 50 CAPSULES - MRP \$268.30 + GST**



# Pink is the New Cure

PinkWave™ was distinctly designed to enhance the efficacy of curing procedures. PinkWave™ is equipped with Vista Apex's patented QuadWave™ Technology which leverages four different wavelengths to ensure peak performance. Along with the standard blue light, PinkWave™ also employs red, near infrared (NIR) and UV light, which make up the only pink curing light on the market. This innovative approach to curing has never been achieved before and results in both doctor and patient benefits.



### Minimised Shrinkage

Powered by Quadwave technology



### Built in Transilluminator

Aids in detection of cracks, caries, interproximal decay



### Increased Depth of Cure

8.5mm



### Distance Forgiveness

Designed to deliver the power needed, regardless the distance from the restoration



### Collimated Beam

Ensures focus from various distances to increase depth of cure



#### V-90734 Pinkwave Curing Light Complete Kit

- 1 x Cordless Curing Light
- 1 x Induction Charging Base
- 1 x Power Cord
- 5 x Autoclavable Light Shields
- 100 x Barrier Sleeves
- 3 x Pair Protective Eyewear

Normally \$3,210.00

**NOW \$2,889.00<sup>ea</sup>**

#### Disposable Barrier Sleeves (100)

**V-90752-2 \$54.90**  
Usually \$61.00



### Articulating paper



This double sided articulating paper is useful to check the occlusion and identify high points in a restoration. Each paper is double sided and 40 microns thick for convenience of use. Packed as 12 sheet booklets, Box of 10 booklets.

- Available in 2 colours: Blue & Red


Normally \$26.50 **\$23.80** ea

### Straight Bite paper holder



**ME-4583** 155mm  
Normally \$40.80 **\$36.70** ea

### Articulating Paper Forceps Miller



**PT-T046**  
Normally \$101.00 **\$90.00** ea

## CONTAC EZ<sup>®</sup> Ultimate Interproximal Solution

### Subgingival Strip System

This system is specially designed to access the subgingival area to make overhang removal simple while preventing damage to patients' gingival tissue and lips. In addition, these ultra-narrow strips are perfect for cleaning out cement debris in the subgingival area of the interproximal space after crown cementation.

#### Cut



**DA-33708**  
Blue Narrow Strip

An ultra-thin 0.07mm double sided serrated strip, without diamonds making it easy to reach and remove subgingival overhangs without any risk of affecting the contact area.

#### Trim



**DA-33408**  
Gold Narrow Strip

A 0.10 mm thick single sided medium fine diamond strip to trim and contour a composite filling.

#### Polish



**DA-33608**  
Copper Narrow Strip

An ultra-thin (0.06 mm) ultra-fine single sided diamond strip to clean and polish a subgingival composite restoration.

**8**  
Pieces

**\$149.40**  
Usually \$166.00

### Subgingival Kit

2x Blue, 3x Gold, 3x Copper

**\$149.40**  
Usually \$166.00 **DA-33108**



**8**  
Pieces




FINISHING & POLISHING SYSTEMS

# OneGloss

One-step 2 in 1 finisher and polisher for composite & enamel.


When used at slow speed with water spray; a firmer pressure will remove composite while softer pressure will polish to a satin finish.

### Cup




Occlusal, Cusp, Incisal

### Midi




Labial, Lingual, Buccal, Occlusal

### IC



Labial, Lingual, Buccal, Interproximal

### Mini



Fissure, Cervical



Item	Code	Usually	
One Gloss Mandrels 6 Pack	SH-0435	\$38.80	<del>\$34.92</del>
Midi Silicone 25 Pack	SH-0171	\$42.80	<del>\$38.52</del>
CA Mini Point Refill 12 Pack	SH-0186	\$42.80	<del>\$38.52</del>
CA Cup Refill 50 Pack	SH-0181	\$85.70	<del>\$77.13</del>
CA Midi Point Refill 50 Pack	SH-0182	\$85.70	<del>\$77.13</del>
CA IC Disc Refill 50 Pack	SH-0183	\$85.70	<del>\$77.13</del>
CA Set 20 of ea + 3 Mandrels	SH-0180	\$102.00	<del>\$91.80</del>



### CA One Gloss Mandrel


For One Gloss Polishers

# Super-Snap

Finishing and polishing disc system.


Designed for easy, fast and safe contouring, finishing and polishing of all micro-filled and hybrid composites.

### Coarse



Contouring

### Medium




Finishing

### Fine



Polishing

### Super Fine



Super Polishing

Safe Side Down	Safe Side Down (m)
SH-L506	SH-L507

Super-Snap Refills 50x  
Diameter 12mm OR Mini 8mm (m)

**\$38.50**

Usually \$42.80

Safe Side Down	Safe Side Up	Double Sided	Double Sided (m)
SH-L508	SH-L509	SH-L528	SH-L522

Double Sided	Double Sided (m)	Double Sided	Double Sided (m)
SH-L511	SH-L531	SH-L512	SH-L532

### Xtreme Mini Kit CA

48 x standard disks, Composite Fine Midi-Point CA, CA Mandrels

**\$57.80**

Usually \$64.30

SH-0509




# SuperBuff

Paste-Impregnated felt disks for direct restorations.



Able to super polish multiple surfaces in a single bound without the mess of polishing pastes.

### Paste-Impregnated



1. Wet the surface to be polished & moisten the SuperBuff disk with water

2. Apply uniform and feather light pressure for polishing



Assorted Set 25x Discs, 2x Mandrels

SH-0535

**\$53.20**


Usually \$59.20

Refills 50x Discs

SH-0536

**\$72.50**

Usually \$80.60



### CA Super-Snap Mandrel

For Super-Snap & Super Buff Discs

Super-Snap Mandrels CA x6

SH-0439

**\$34.90**

Usually \$38.80

# Contours™

Coarse Anatomy Trimmers

- 90 micron aluminium oxide coarse anatomy trimmers
- Trim composite quickly, without significant degradation
- Autoclavable
- Available in 4 shapes perfectly suited for every anterior or posterior restoration

Assorted Kit

**CC-078912**

12 Pack

**\$61.20**

Usually \$68.00



- #01 CC-078201
- #04 CC-078204
- #05 CC-078205
- #06 CC-078206
- #08 CC-078208

**Clinician's Choice**  
Designed for better dentistry.



20 Packs  
**\$101.70**  
Usually \$113.00

# A.S.A.P.®

Diamond Polishing Brushes

**Clinician's Choice**  
Designed for better dentistry.

# A.S.A.P.®

All Surface Access Polishers



A.S.A.P.® Diamond Polishing Brushes quickly produce a high gloss polish on both incline planes, ridges, pits, and grooves, all with one size and shape. Their diamond impregnated bristles spread over cusp tips and marginal ridges, while the individual bristles polish the narrowest of grooves without altering anatomy.

**Trial Kit**

- 1x Pre-Polisher
- 1x Final-Polisher

**CC-133159**

**\$63.00**

Usually \$70.00

**3 Pack Pre-Polishers**

**CC-133156**

**\$98.00**

Usually \$109.00

**3 Pack Final-Polishers**

**CC-133157**

**6 Pack Polishers**

- 3x Pre-Polishers
- 3x Final-Polishers

**CC-133158**

**\$188.00**

Usually \$209.00


**A.S.A.P. Diamond Brush & Wheel System Kit**

- 2x Pre-Polisher Wheels (Large)
- 2x Pre-Polisher Wheels (Small)
- 2x Final High Shine Wheels (Large)
- 2x Final High Shine Wheels (Small)
- 2x Pre-Polisher Diamond Brushes
- 2x Final Polisher Diamond Brushes

**CC-133155**


**\$327.60**

Usually \$414.00



**PRE**  
Polish

**FINAL**  
Polish



**2 Pack Small**

- 1x Pre-Polisher
- 1x Final-Polisher

**CC-130910**

**\$58.00**

Usually \$64.00

**2 Pack Large**

- 1x Pre-Polisher
- 1x Final-Polisher

**CC-130911**

**\$73.00**

Usually \$81.00

**3 Pack Small Pre-Polishers**

**CC-131011**

**\$83.00**

**3 Pack Small Final-Polishers**

**CC-132013**

Usually \$92.00

**3 Pack Large Pre-Polishers**

**CC-131012**

**3 Pack Large Final-Polishers**

**CC-132014**

**\$102.60**

Usually \$114.00



**centrix**  
MAKING DENTISTRY EASIER

# FluoroDose®

5% Sodium Fluoride Varnish with Xylitol



## Dip. Mix. Brush. Throw away.

FluoroDose contains 5% sodium fluoride in a convenient single-use Lollipop®. Each unit contains one single dose (0.3mL) of FluoroDose plus one Benda® Brush applicator. Simply peel back the cover, dip in the brush, mix and apply.

## Description

- Contains Xylitol
- Unit-dose delivery allows for ease of use and clean-up
- Swipe-on procedure is easier and safer to use than gels and rinses, done in less than a minute
- Dries in seconds upon contact with saliva
- Remains on tooth for 6 to 8 hours for optimum fluoride uptake
- Available in 5 great-tasting flavours: Caramel, Bubble Gum, Mint, Cherry and Melon
- 5% sodium fluoride (NaF) or 22,600 ppm of fluoride
- Benda Brush applicators also allow you to mix the varnish prior to application, guaranteeing a consistent level of fluoride in every application
- No more messy gels and uncomfortable trays; no more inherent worry of fluoride ingestion

## Allergen Risk Information

- FluoroDose is lactose free, gluten free and nut free
- Product is free of all dairy, egg, gluten, soy, peanut, tree nuts, tree nut derivatives, fish/shellfish and seed products
- FluoroDose contains rosin, beeswax and red dye 40

## 120 PACK

- |                             |                          |
|-----------------------------|--------------------------|
| <b>CE-360086</b> MINT       | <b>CE-360110</b> CHERRY  |
| <b>CE-360087</b> BUBBLE GUM | <b>CE-360137</b> CARAMEL |
| <b>CE-360107</b> MELON      |                          |
- Normally \$326.40 **\$293.70** ea

## 5 FLAVOUR MIXED PACK

- CE-700001** FLUORIDE VARNISH 40 PIECE
- Normally \$112.20 **\$100.90** ea

Perfect for private practice



## Flavours

- Cherry 
 Mint 
 Caramel 
 BubbleGum 
 Melon



# DIRECTA

# PractiPal®

Quiet, light-weight and colourful, PractiPal is a pleasure to work with! It is made of highly durable plastic that does not absorb heat in steam autoclaves and washer disinfectors, aiding the cleaning process and saving energy costs. PractiPal Tray System includes trays with bur stands, file stands, consumables

and accessories that can easily be secured onto the trays throughout the entire cleaning and disinfection process. Both the tray clamp and bur stand/file holder are easy to open. Once prepared, a tray configuration can be cycled over and over again.



## PractiPal® Mini Trays

Pink DA-1150-56  
Grey DA-1150-50



PractiPal in a smaller size holding **5 instruments** - ideal for checkups, for dental hygienists and to stock spare instruments. Dimensions: 80 x 188 x 37mm.

**\$107.10**  
Usually \$119.00

## PractiPal® Half Trays

Lilac DA-1151-06  
Blue DA-1151-13



PractiPal Half Tray suits all set-ups and contains one half tray and one instrument clamp, holding up to **9 instruments**. Dimensions: 140 x 188 x 37mm.

**\$127.80**  
Usually \$142.00


All PractiPal trays available in the following colours on request: **Black/**



## PractiPal® Full Tray Complete Sets



Mint Green  
DA-1151-97



PractiPal Full Tray with place for **9 instruments**  
Complete Set contains: Tray (2 pcs), clips (4 pcs), Instrument Clamp (1 pc), Bur/Endo stand for 16 burs or 28 files (1 pc), MultiBlock (1 pc), MultiDappen (5 pcs), MultiPoint (5 pcs), Waste Cup (5 pcs), Foam Large (1 pc) and Foam Small (5 pcs).  
Dimensions: 280 x 188 x 37mm.

**\$453.60**  
Usually \$504.00

See page 10 for PractiPal® Bur and File Stands



## Fliptop™ Cassettes

PDT's FlipTop™ Cassettes are made with the highest quality, passivated, surgical grade stainless steel and a unique Teflon double hinge that flips up, over, and under for a much smaller footprint. Save time and prolong the life of your instruments while maintaining maximum efficiency and organisation.

PDT cassettes are available in configurations that will hold from 5 to 22 instruments.

**Double Racks (D)**  
**Fixed (F)**  
**Liftout (L)**  
**Utility (U)**

**Contact your HeathCare Essentials Representative** for more information on the PDT Cassette range.



PT-T005A-F

### Fliptop Cassette


A Fixed (F) Rack  
5 Instrument

**\$172.00**

Usually \$192.00

### Fliptop Cassette

B Liftout (L) Rack 7 Instrument



PT-T007B-L

**\$241.00**

Usually \$268.00

### Fliptop Cassette

C Fixed (F) Rack 9 Instrument




PT-T009C-F

**\$284.00**

Usually \$316.00

### Fliptop Cassette

D Double (D) Rack 18 Instrument



PT-T009D-D

**\$315.00**

Usually \$350.00

### Fliptop Cassette

D Liftout (L) Rack 9 Instrument



PT-T009D-L

**\$324.00**

Usually \$360.00





## Organise Your Time.

Choose from 5 beautiful silicone insert colours, specifically developed for medical use. Safe non-toxic and durable. The special shape allows an extremely safe handling of surgical instruments during the washing and sterilisation phases.



The exclusive Gammatrix locking system is practical and easy to use. It keeps the tray securely closed in all conditions. The open corners facilitate the outflow of water and residue reducing drying times.




- A. Open corners allow residuals and water to flow out from the tray reducing drying time
- B. Specifically developed for medical purposes, this atoxic silicon doesn't release dangerous agents on instruments. It respects the CE Food's normative
- C. Safe and easy joint system, totally hygienic design
- D. Open spaces assure the perfect washing process and safe sterilisation.
- E. V-form upper silicone provides secure storage for instruments of differing sizes




**\$90.90**  
Usually \$101.00

- Blue Silicone**  
ME-978-BL
- Red Silicone**  
ME-978-RO
- Yellow Silicone**  
ME-978-GI
- Green Silicone**  
ME-978-VE
- Black Silicone**  
ME-978-NE




**\$179.00**  
Usually \$198.90

- Blue Silicone**  
ME-981-19-BL
- Red Silicone**  
ME-981-19-RO
- Yellow Silicone**  
ME-981-19-GI
- Green Silicone**  
ME-981-19-VE
- Black Silicone**  
ME-981-19-NE




**\$137.70**  
Usually \$153.00

- Blue Silicone**  
ME-980-9-BL
- Red Silicone**  
ME-980-9-RO
- Yellow Silicone**  
ME-980-9-GI
- Green Silicone**  
ME-980-9-VE
- Black Silicone**  
ME-980-9-NE



**\$206.50**  
Usually \$229.50

- Blue Silicone**  
ME-981-20-BL
- Red Silicone**  
ME-981-20-RO
- Yellow Silicone**  
ME-981-20-GI
- Green Silicone**  
ME-981-20-VE
- Black Silicone**  
ME-981-20-NE



**\$137.70**  
Usually \$153.00

- Blue Silicone**  
ME-980-10-BL
- Red Silicone**  
ME-980-10-RO
- Yellow Silicone**  
ME-980-10-GI
- Green Silicone**  
ME-980-10-VE
- Black Silicone**  
ME-980-10-NE

Scan to learn more about  
**GAMMAFIX**  
by MEDESYS



# ANATOMIC TOOTH FORCEPS

## UPPER MOLARS



ME-2500-17



ME-2500-17

Upper Molars, Right



ME-2500-18

Upper Molars, Left



ME-2500-18-A

Upper Molars Either Side



ME-2500-89

Upper Molars, Right



ME-2500-90

Upper Molars, Left



ME-2500-67

Upper Third Molars and  
Upper Wisdoms Either Side



ME-2500-67-A

Upper Third Molars and  
Upper Wisdoms Either Side



ME-2500-19

Upper Third Molars and  
Upper Wisdoms Either Side




ME-2500-54

For Separating Upper Molars



ME-2500-55

For Separating Upper  
Molars




ME-2500-94

Upper Molars Right



ME-2500-95

Upper Molars Left



ME-2500-96

Upper Molars Right



ME-2500-97

Upper Molars Left

ME-2500-96 and  
ME-2500-97 have  
slightly shorter points  
than ME-2500-94 and  
ME2500-95.

Normally \$153.00  
**\$137.70**  
ea

\* Some items may be on backorder.

# LOWER MOLARS



**ME-2500-73**

(1:1 scale)



**ME-2500-73**

Lower Molars




**ME-2500-73-S**

Lower Molars



**ME-2500-73-L**

Lower Molars



**ME-2500-21**

Lower Molars




**ME-2500-22**

Lower Molars




**ME-2500-148**

Balding,  
Lower Molars



**ME-2500-23**

Lower Molars, Right




**ME-2500-24**

Lower Molars, Left



**ME-2500-20**

Lower Third Molars



**ME-2500-79**

Lower Third Molars and  
Lower Wisdoms Either Side



**ME-2500-79-A**

Lower Third Molars

Normally \$153.00  
**\$137.70** ea

# COW HORN BEAK



**ME-2500-86**

Lower Molars Decayed  
or Broken Down Crowns



**ME-2500-86C**

Lower Molars Decayed  
or Broken Down Crowns



**ME-2500-91**

Rowney, For Separating  
Lower Molars



**ME-2500-87**

Lower Molars Decayed  
or Broken Down Crowns

Normally \$153.00  
**\$137.70** ea

\* Some items may be on backorder.

# Every Indication. Every Workflow.

Ideal tray material/light body synergy to ensure light body is not overdisplaced, resulting in ultimate impression accuracy.

Branched resin 3-dimensional cross linking results in an ultra high tear strength.

Hydroactive chemistry gives consistent results in all oral conditions.

9 Technique designed viscosities.

Doctor-friendly long working time  
Patient-friendly short intraoral cure time.

## Affinity™



Hydroactive Impression Material

With a unique branched resin, more refined filler and an even distribution of surfactant, Affinity™ Hydroactive Impression Material is clinically proven to provide consistently accurate impressions.



### Spectrum Packs

Spectrum Pack supplied in quantities of 12, 24, or 36.

You can choose any combination of Affinity Singles to make up a spectrum pack.

12 singles	\$528.00
24 single	\$990.00
36 singles	\$1425.60










For codes: see SS under Pack Options on table opp.:



# Affinity™

Hydroactive Impression Material

## Technique Designed Viscosity Guide

Affinity Viscosity		INTENDED USE	WORKING TIME	MINIMUM INTRAORAL SET TIME	TRAY OPTIONS	PACK OPTIONS
LIGHT BODY (HIGH FLOW) REG. SET		Flows on its own under the hydraulic pressure of the tray material and flows effortlessly into the sulcus. Our most popular light body impression material.	1:45 min.	2:30 min.*	<ul style="list-style-type: none"> <li>• Quad-Trays®</li> <li>• HeatWave®</li> <li>• BorderLock®</li> </ul>	2 x 50ml CC-024420 \$139.00 SS CC-019510
LIGHT BODY (HIGH FLOW) FAST SET			1:10 min.	1:45 min.**		2 x 50ml CC-024425 \$139.00 SS CC-019525
LIGHT BODY (REGULAR FLOW) REG. SET		Thixotropic light body impression material, designed for the "air assisted" placement technique.	1:45 min.	2:30 min.*	<ul style="list-style-type: none"> <li>• Quad-Trays®</li> <li>• HeatWave®</li> <li>• BorderLock®</li> </ul>	2 x 50ml CC-024410 \$139.00 SS CC-019510
LIGHT BODY (REGULAR FLOW) FAST SET			1:10 min.	1:45 min.**		2 x 50ml CC-024415 \$139.00 SS CC-019515
LIGHT BODY XL (XTRA LIGHT)		Most fluid viscosity light body. Recommended for use with healthy soft tissue when it is difficult to impress margins.	1:45 min.	2:30 min.*	<ul style="list-style-type: none"> <li>• Quad-Trays®</li> <li>• HeatWave®</li> <li>• BorderLock®</li> </ul>	2 x 50ml CC-024426 \$139.00 SS CC-019526
HEAVY BODY REG. SET		Ideal for both quadrant and full-arch impressions and ensures a perfect blend with the light body, preventing overdisplacement.	1:45 min.	2:30 min.*	<ul style="list-style-type: none"> <li>• HeatWave®</li> <li>• BorderLock®</li> </ul>	2 x 50ml CC-023400 \$139.00 SS CC-013912
HEAVY BODY FAST SET			1:10 min.	1:45 min.**		2 x 50ml CC-023405 \$139.00 SS CC-013952
INFLEX MAXIMUM SUPPORT TRAY MATERIAL REG. SET		Specifically designed for use with a dual-arch impression technique. Rigid tray material designed to prevent impression distortion.	1:45 min.	2:30 min.*	<ul style="list-style-type: none"> <li>• Quad-Trays®</li> <li>• HeatWave®</li> <li>• BorderLock®</li> </ul>	2 x 50ml CC-022400 \$139.00 SS CC-012912
INFLEX MAXIMUM SUPPORT TRAY MATERIAL FAST SET			1:10 min.	1:45 min.**		2 x 50ml CC-022405 \$139.00 SS CC012952
QUICK BITE BITE REGISTRATION MATERIAL		Mousse-like bite registration with an ideal flow to capture accurate detail. Carvable, ultra-high 90 durometer cure eliminates distortion during prosthetic fabrication.	0:15 min.	0:45 min.		2 x 50ml CC-026400 \$139.00 SS CC-016912
MONOPHASE		Single phase impression material designed for implant impressions, complete/partial denture impressions, and fabrication of clear aligners. Monophase features a gum-like consistency, high durometer cure and high tear strength while also being sterilisable and dimensionally stable.	2:15 min.	2:45 min.	<ul style="list-style-type: none"> <li>• HeatWave®</li> <li>• BorderLock®</li> </ul>	2 x 50ml CC-025400 \$139.00 SS CC-019540
AFFINITY CRYSTAL CLEAR MATRIX IMPRESSION MATERIAL		Water-clear silicone matrix material designed for easy fabrication of temporary veneers. Flowing effortlessly into embrasure areas, Crystal captures extreme detail. Its transparency allows control of flash formation, seating, and easy detection of voids and air bubbles.	0:45 min.	1:30 min.	<ul style="list-style-type: none"> <li>• Clear Impression Trays</li> </ul>	2 x 50ml CC-021400 \$190.00 6 x 50ml CC-019916 \$416.00 SS CC-010912
AFFINITY PUTTY		Non-greasy and satiny smooth putty designed to ensure a feathered edge at the interface with wash material.	1:45 min.***	3:15 min.	<ul style="list-style-type: none"> <li>• HeatWave®</li> <li>• BorderLock®</li> </ul>	CC-027200 \$237.00

\*2:30 min. reflects the minimum intraoral set time. Our recommendation is to allow 3:00 min. for a full intraoral set.

\*\*1:45 min. reflects the minimum intraoral set time. Our recommendation is 2:00 min. for full intraoral set.

\*\*\*Available working time for Putty is in addition to a 0:30-0:45 second mix time.

These procedure times reflect ideal conditions. Fluctuations in temperature and humidity may reduce or extend these times.

## Luxatemp Fluorescence

The temporary crown and bridge material for natural aesthetics.



- Natural fluorescence in six shades
- Low polymerisation shrinkage
- High fracture and abrasion resistance
- Outstanding colour stability
- Low allergy potential
- Suited for long-term temporaries



1x 76g Cartridge  
15x Automix-Tips

**\$216.60**  
Usually \$240.70

**A1** DM-110585  
**A2** DM-110586  
**A3** DM-110710

**A3.5** DM-110587  
**B1** DM-110588  
**BL** DM-110589

## Luxatemp Automix Tips

Mixing tips for Luxatemp Automix Cartridges


- Mixing Ratio 1:1
- Blue / Orange

45 Pack

DM-110409

**\$64.20**

Usually \$71.40



## Luxatemp Glaze & Bond

Surface sealer and repair material for temporary crowns and bridges.

- High gloss shine without polishing
- Excellent bond with repair materials
- Extraordinary abrasion resistance
- Low surface roughness



**\$72.50**

Usually \$80.60

DM-211927

5 ml Bottle

## TempoCem NE luting cement

Temporary luting of crowns, bridges, inlays and onlays.

TempoCem NE is your self-mixing, classic, temporary cement without eugenol. It is ideally suited to temporary luting of temporary and permanent crowns, bridges, inlays and onlays, not least due to its very thin film thickness.



DM-212104

**\$101.00**  
Usually \$112.20

**Smartmix**

2x Smartmix Syringes @ 5 ml  
20x Smartmix Tips Short

## Tray Adhesive For A-Silicones

Tray adhesive for all VPS materials

DMG Tray Adhesive for A-Silicones is formulated to offer optimal adhesive between addition-silicone-based impression materials (VPS) and minimally retentive or non retentive metal, plastic or resin trays.

It is ideal for use with all impression materials, especially Honigum. Available in convenient 10 ml bottle.



- Promotes the bond of VPS-based impression materials in the impression tray
- Easy to apply and handle
- Short setting time after application
- Blue colour visual control

10 ml Bottle

DM-909394

**\$45.90**

Usually \$51.00

**DMG**  
Dental Milestones Guaranteed



**GERMIPHENE®**  
**Tres Clean**  
**Tray Cleaner**

A mild cleaning formula that leaves trays and instruments sparkling clean. Safe for all plastics and non-oxidising metals. Non-rusting 1 litre bottle of concentrate replaces 20 litre of ready-to-use solution.

**\$117.50**  
Usually \$130.60

GP-HTR 1 Litre Bottle



**PULPDENT®**  
**Tray Adhesive**

For use with all impression materials. Place a thin application on the impression tray.

**\$32.40**  
Usually \$36.00

PL-DDTA2 30ml Bottle



**Alginmax Alginate**

High precision alginate for impressions with chromatic phase indicator.

Thanks to the MAP technology (modified atmosphere packaging) the product shelf life is 6 years.

The excellent dimensional stability allows casting of the impressions after 5 days (if correctly stored in sealed container).

**SETTING TIMES (at 23°C):**  
Mixing .....35sec  
Working.....1min 20sec  
Min Setting in mouth.... 40sec  
Initial Setting.....1min 20sec  
Total Setting ..... 2min

Buy

12

**20%**  
OFF

Buy

24

**30%**  
OFF



**MEDESÝ** ITALY **Impression Tray Kits**

**Perforated with Rim**



12 pcs (6 upper, 6 lower)  
XXS - XL Upper  
XXS - XL Lower

**ME-6000-KIT**  
**\$179.00**  
Usually \$198.90

**Solid with Rim**



12 pcs (6 upper, 6 lower)  
XXS - XL Upper  
XXS - XL Lower

**ME-6005-KIT**  
**\$179.00**  
Usually \$198.90

**Edentulous Perforated Special Shape**



11 pcs (5 upper, 6 lower)  
XXS - L Upper  
XXS - XL Lower

**ME-6011-KIT**  
**\$183.60**  
Usually \$204.00

**Edentulous Solid with Rim**



11 pcs (5 upper, 6 lower)  
XXS - L Upper  
XXS - XL Lower

**ME-6009-KIT**  
**\$183.60**  
Usually \$204.00

## Vista Green Impression Trays

- Superior visibility and greater accuracy
- 9 sizes for easy anatomy matching



<b>V-311005</b> Small Upper		
<b>V-311006</b> Small Lower		
<b>V-311003</b> Medium Upper		
<b>V-311004</b> Medium Lower		
<b>V-311001</b> Large Upper	<b>V-311008</b> UL/LR	<b>V-311009</b> Anterior
<b>V-311002</b> Large Lower	12 pack	
<b>V-311007</b> UR/LL	Normally \$29.00	<b>\$26.00</b> ea

## Clear Impression Trays

Clear, non-perforated impression trays. Designed for use with Template clear to create a clear matrix for veneer or crown tempor.



<b>Full arch upper (medium)</b>	<b>Universal Anterior (upper/lower)</b>
<b>CC-036601</b> 6 pack	<b>CC-036603</b> 6 pack
<b>CC-036251</b> 25 pack	<b>CC-036253</b> 25 pack
<b>Full arch lower (medium)</b>	
<b>CC-036602</b> 6 pack	6 pack
<b>CC-036252</b> 25 pack	Normally \$25.00
	25 pack
	Normally \$69.00
	<b>\$22.50</b> ea
	<b>\$62.10</b> ea


## Unident Trays

- Extra rigid plastic
- Distortion free
- Designed for maximum retention
- Pack of 12



<b>GB-IMP001</b> Small Upper	
<b>GB-IMP002</b> Small Lower	
<b>GB-IMP003</b> Medium Upper	
<b>GB-IMP004</b> Medium Lower	
<b>GB-IMP005</b> Large Upper	
<b>GB-IMP006</b> Large Lower	
	12 pack
	Normally \$28.60
	<b>\$25.70</b> ea

## Edentulous Trays



<b>GB-U3ESU</b> Small Upper	
<b>GB-U3ESL</b> Small Lower	
<b>GB-U3EMU</b> Medium Upper	
<b>GB-U3EML</b> Medium Lower	
<b>GB-U3ELU</b> Large Upper	
<b>GB-U3ELL</b> Large Lower	
	12 pack
	Normally \$28.60
	<b>\$25.70</b> ea

## Heatwave<sup>®</sup>

Customisable Impression Trays

With HeatWave, the clinician can easily self-customise the thermoplastic trays by hand after immersing them in water for 60 seconds at 70°C. The tray can then be shaped by hand and moulded to a patient's particular arch size and contour. After adaption, the tray will assume its normal rigidity at room temperature or by running under cold water, allowing for a precise impression.



**Clinician's Choice**  
Designed for better dentistry.

Size 42	Upper Tray	CC-188142
Size 43	Upper Tray	CC-188143
Size 44	Upper Tray	CC-188144
Size 45	Upper Tray	CC-188145
Size 51	Lower Tray	CC-188151
Size 52	Lower Tray	CC-188152
Size 53	Lower Tray	CC-188153
Size 54	Lower Tray	CC-188154

Intro Kit  
x4 of each size, Measuring Calipers & Sizing Chart

**\$115.20**  
Usually \$128.00 **CC-188908**




10 Packs  
**\$80.00**  
Usually \$89.00



## QUAD-TRAY Xtreme



- Wide, customisable arch design to fit any arch shape or size
- Thin distal bar, to prevent retro molar pad impingement
- Large handle for easy handling
- Ultra low to prevent palatal impingement



x50 Posterior  
 Normally \$249.00  
**\$224.00** ea  
 CC-048050

## TempTrays™


Designed for accuracy, TempTray™ is a distortion-free, disposable, metal temporary impression tray with extreme versatility. Available in anterior and posterior variations for use in either upper or lower arch. The sidewalls of the tray are high to encompass virtually any crown length.



x50 Anterior Normally \$156.00  
**\$140.40** ea CC-046150  
 x50 Posterior Normally \$121.00  
**\$108.90** ea CC-046050

## QUAD-TRAY XTRA LONG XL

- High side walls to accommodate all tooth lengths
- Perforated walls to retain impression material
- Customisable arch to fit all patients
- Wide handle for easy mouth removal and handling
- Recyclable aluminium



CC-041400  
 40 per pk - XL

Normally \$254.00  
**\$228.60** ea

## 3.n.One™ Impression Trays

For Dual-Arch Impression Technique

- Unique undercut ridges lock impression materials in place
- High rigidity with non-tear, flexible webbing ensures maximum accuracy
- Five different shapes to meet all applications



CE-360011 Anterior 35 pk  
 CE-360012 Posterior 50 pk  
 CE-360013 Quad 35 pk

CE-360014 Mixed 40 pk  
 Normally \$203.00  
**\$182.70** ea

## Redi-Lock Impression Trays

While stocks last

Redi-Lock Impression Trays are made of advanced medical-grade polystyrene, which is designed to be disposable. They will not bend or flex and always provide sharp detail. Can be easily customised with hot water, and feature sturdy handles for better control. Sold individually or in bags of 50.



MO-2019601 Child Small Lower Red  
 MO-2019602 Child Small Upper Red  
 MO-2019603 Child Medium Lower Orange  
 MO-2019604 Child Medium Upper Orange

MO-2019605 Child Large Lower Blue  
 MO-2019608 Adult Small Upper Green  
 MO-2019609 Adult Medium Lower Yellow  
 MO-2019610 Adult Medium Upper Yellow

MO-2019611 Adult Large Lower Purple

Normally \$75.00  
**\$67.50** ea



# Tips on Tips

## Universal Tips

## Tissue Management Tips

## Endodontic Tips

## Periodontal Tips

### Your Total Tip Solutions

Dental syringes are some of the most common pieces of equipment used by dental professionals in a wide range of dental procedures, from the most routine to the most complicated surgical procedures. At Vista Apex, we offer complete syringes, dental syringe tips (including bendable tips, dripless tips, brush tips, and more), and accessories like tip covers.

With many of these products, we offer the option to buy in bulk so your dental practice is always stocked with dental syringes and syringe tips. As with all our dental products, we aim to deliver the highest quality materials to make sure your practice is treating patients right.

**Tissue Management**

**Stat-Flo™ Brush Tip**

- Padded-end brush tip



**20 cup**  
V-312102  
**100 cup**  
V-312100

**Designed for:**  
Hemostatic agents

**MST™ Universal Tip**

- Material saving technology®



**100 cup**  
22ga  
V-315310

**Designed for:**  
Flowable composites, pulp capping, etchants, & sealants

**Spira-Flo™ Brush Tip**

- Nylon brush fibers for controlled delivery
- Internal spiral
- Chamber minimizes bubbles



**100 cup**  
V-321400

**Designed for:**  
Etchants and resins

**Sol-u-Flo™ Applicator Tip**

- 19ga tip with adjustable brush fiber length
- Curved cannula for optimal entry



**20 cup**  
V-312202  
**100 cup**  
V-312200

**Exacta-Flo™ Brush Tip**

- Adjustable, soft brush fibers minimize bubbles



**100 cup**  
V-321610

**Designed for:**  
Aqueous solutions such as Caries Indicator, Silane, and Chlorhexidine

**Micro-Flo™ Tip**

- Polyimide tip for microscope enhanced dentistry



**50 cup**  
26ga  
V-406626  
27.5ga  
V-406627

**Designed for:**  
Flowable composites, pulp capping, etchants, & sealants

**Volu-Flo™ Tip**

- Delivers large volumes of material (18ga)



**100 cup**  
V-321100

**Designed for:**  
Viscous materials, pastes & putties

**1/2" Bendable Applicator Tip**

- Blunt end, dead-soft needle resists crimping



22ga  
**20 cup**  
V-315202  
**100 cup**  
V-315222

25ga  
**20 cup**  
V-315205  
**100 cup**  
V-315225

**Designed for:**  
Etchants, resins, and flowable composites

**Blue-Flo™ Tip**

23ga soft plastic tip



**20 cup**  
V-316202  
**100 bag**  
V-316222

**Designed for:**  
Irrigation, aeration & impression materials


SCAN HERE  
For More Information



VISTA APEX  
vistaapex.com  
877.418.4782

**Appli-Vac™ 3/4" Bendable Tip**

- Blunt end, dead-soft needle resists crimping



**100 cup**  
16ga  
V-315116

**100 cup**  
18ga  
V-315118

**100 cup**  
19ga  
V-315119

**100 cup**  
20ga  
V-315120


**100 cup**  
22ga  
V-315122

**100 cup**  
25ga  
V-315125

**Designed for:**  
Sealants, flowable composites, etchants, resins, primers

**Pre-Bent Tip**

- Blunt end pre-bent needle
- Opaque photo sensitive hub (22/25ga)



**100 cup**  
16ga  
V-315316

**100 cup**  
20ga  
V-315320

**100 cup**  
22ga  
V-315322

**100 cup**  
25ga  
V-315325

**Designed for:**  
Etchants, resins and flowable composite



**MST Visco-Tip™**

- Bendable, flexible 25ga tip




**20 cup**  
V-312302

**50 cup**  
V-312300

**Designed for:**  
Viscous materials such as bioceramic sealers and calcium hydroxide

**Voyager Irrigating Tip**

- Dual side ports for safety and efficacy
- Enhanced flexibility to navigate canals
- Silicone stopper for depth reference




<b>17mm</b> 27ga V-315727	<b>17mm</b> 29ga V-315729	<b>17mm</b> 30ga V-315730	<b>17mm</b> 31ga V-315731
<b>27mm</b> 27ga V-315827	<b>27mm</b> 29ga V-315829	<b>27mm</b> 30ga V-315830	<b>27mm</b> 31ga V-315831

**Designed for:**  
Endo irrigating solutions

**Appli-Vac™ 1" Bendable Irrigating Tips**

- Slotted and side-vented



<b>20 cup</b> 23ga V-315023	<b>20 cup</b> 27ga V-315027	<b>20 cup</b> 30ga V-315030
<b>100 bag</b> 23ga V-315523	<b>100 bag</b> 27ga V-315527	<b>100 bag</b> 30ga V-315530

**Designed for:**  
Endo Irrigating solutions (27 ga / 30 ga)  
Perio Irrigation (23 ga)

**Secure-Lock™ Appli-Vac™ 1" Irrigating Tip**

- Slotted and side-vented
- Secure-Lock™ threads lock tip in place




<b>80 bag</b> 23ga V-314523	<b>80 bag</b> 27ga V-314527	<b>80 bag</b> 30ga V-314530
-----------------------------------	-----------------------------------	-----------------------------------

**Designed for:**  
Endo Irrigating solutions (27 ga / 30 ga)  
Perio Irrigation (23 ga)

**Vista-Probe™ Irrigating Tip**

- Closed-end for side-port delivery
- Lateral delivery of irrigants reduces apical pressure




<b>20 cup</b> 23ga V-315123	<b>20 cup</b> 27ga V-315127	<b>20 cup</b> 30ga V-315130
<b>100 bag</b> 23ga V-315623	<b>100 bag</b> 27ga V-315627	<b>100 bag</b> 30ga V-315630

**Designed for:**  
Endo Irrigating solutions (27 ga / 30 ga)  
Perio Irrigation (23 ga)

**Secure Lock™ Vista-Probe™ Tip**

- Closed-end for side-port delivery
- Lateral delivery of irrigants reduces apical pressure
- Secure-Lock™ threads lock tip in place



<b>80 bag</b> 23ga V-314623	<b>80 bag</b> 27ga V-314627	<b>80 bag</b> 30ga V-314630
-----------------------------------	-----------------------------------	-----------------------------------

**Designed for:**  
Endo Irrigating solutions (27 ga / 30 ga)  
Perio Irrigation (23 ga)

**Elasti-Vac™ Tip**

- Narrow, flexible tip easily navigates curved canals



**50 cup**  
V-406150

**Designed for:**  
Endodontic evacuation

**Versa-Brush™**

- Rotational utility brush



**40 cup**  
V-204300

**Designed for:**  
Cleaning post holes and coronal access openings

**NiTi SuperFlex™ Tip**

- AUTOCLAVABLE. Slotted and side-vented, 30ga for safe irrigation
- 2.5x more flexible



<b>17mm (short)</b> 6 cup V-405107	<b>25mm (long)</b> 6 cup V-405108	<b>20 cup</b> V-406320
<b>Mixed 6 cup</b> V-405106	<b>50 cup</b> V-406350	

**Designed for:**  
Endo irrigating solutions

**Periodontal Pocket-Probe™**

- Closed-end with multiple perforations for thorough flushing



**100 cup**  
21ga  
V-316102

**Designed for:**  
Flushing of perio pockets or inter-dental spaces

### Automix Dispensing Tips

 <p><b>TIPS BLUE 10:1</b> Suitable for: Temporisation <b>PP-7010</b> 50pk</p>	 <p><b>TIPS GREEN 1:1</b> Suitable for: Impressions <b>PP-7002</b> 50pk</p>
 <p><b>TIPS PINK 1:1</b> Suitable for: Bite and reline material <b>PP-7003</b> 50pk</p>	 <p><b>TIPS YELLOW 1:1</b> Suitable for: Impressions <b>PP-7001</b> 50pk</p>

Normally \$32.60 **\$29.30**<sub>ea</sub>

### Intraoral Tips

### Automix Tips

 <p><b>DM-909201</b> Yellow</p>	 <p><b>DM-909200</b> Yellow Normally \$90.80</p>
 <p><b>DM-909400</b> Clear</p>	<p>• Universal fitting • Available in yellow or clear • 50 per pk</p> <p>Normally \$51.00 <b>\$45.90</b><sub>ea</sub></p>

**\$81.70**<sub>ea</sub>

### Mixstar Mixing Tips

**DM-909567**  
White 50pk

Normally \$155.00 **\$139.50**<sub>ea</sub>

### Smart Mix Tips

Short  
**DM-212040** 50pk  
**\$77.10**<sub>ea</sub>  
Normally \$85.70




Combi  
**DM-212042** 50pk  
**\$154.20**<sub>ea</sub>  
Normally \$171.40

Endo  
**DM-212041** 50pk  
**\$87.20**<sub>ea</sub>  
Normally \$96.90

### Mixing Tips

With Integrated Needle

Indicated for impression material as well as core build-up material.

 <p><b>CC-032252</b> 2.8 Brown 25 Pack</p>	 <p><b>CC-032253</b> 3.2 Blue 25 Pack</p>	 <p><b>CC-032254</b> 4.2 Yellow 25 Pack</p>
--	--	--

Normally \$100.00 **\$90.00**<sub>ea</sub>

### Mixing Tips

Ideal for affinity crystal as well as counter-fit multi-purpose replication silicone.

Indicated for all affinity viscosities, but ideal for tray materials and multi-prep tray material.

 <p><b>CC-033124</b> 7.5 Purple 100 Bag Normally \$205.00</p>	 <p><b>CC-033126</b> 6.5 Aqua, 100 Pack Normally \$226.00</p>
--	--

**\$184.50**<sub>ea</sub> **\$203.40**<sub>ea</sub>

### NeoSEALER Flex Flo Tips

- Integrated stopper for accuracy of placement
- 25.5 gauge flexible tip to reach deep into curved and small canals

<p>20 Pack <b>AV-FFT20</b> <b>\$73.40</b> Usually \$81.60</p>	<p>40 Pack <b>AV-FFT40</b> <b>\$128.50</b> Usually \$142.80</p>
---	---



### NeoSEALER Flo Basic Tips

- Highly flexible to follow canal anatomy without crimping or binding
- 25 gauge flexible and trimmable tip

<p>20 Pack <b>AV-BFT20</b> <b>\$44.10</b> Usually \$49.00</p>	<p>40 Pack <b>AV-BFT40</b> <b>\$67.90</b> Usually \$75.50</p>
---	---



### NeoSEALER Flo Tips

- Bendable tip, without crimping, for better accessibility
- 25 gauge, 25mm cannula length

<p>20 Pack <b>AV-CKT20</b> <b>\$40.40</b> Usually \$44.90</p>	<p>40 Pack <b>AV-CKT40</b> <b>\$60.50</b> Usually \$67.30</p>
---	---



**DIRECTA**

**15% OFF**

**ProphyCare®**

**ProphyCare® Prophy Paste CCS various grits for various needs**

Prophy Paste CCS is second to none in performance, handling and dispensing. The unique non-splattering, homogeneous formula does not dry out, liquify or separate. Thick layers of plaque and stains are easily removed leaving a fresh and smooth minty sensation after each treatment. Packed in easy to dispense, **colour coded tubes** and **single dose disposable cups (144x 2g)**. Available in 4 different grits from coarse to extra-fine. **Contains 0,1% fluoride.**

**"The perfect starter, effective removal of stubborn stains"**



**COARSE Blue RDA 250**

		Usually	NOW
Tube 60ml	DA-7311-13	\$48.00	<b>\$40.80</b>
Paraben Free	DA-7311-16	\$48.00	<b>\$40.80</b>
Single Dose	DA-6901-23	\$215.00	<b>\$182.70</b>

**"Less abrasive, but still effective"**



**MEDIUM Green RDA 170**

		Usually	NOW
Tube 60ml	DA-7311-14	\$48.00	<b>\$40.80</b>
Single Dose	DA-6901-22	\$215.00	<b>\$182.70</b>

**"Gentle cleaning, tissue friendly stain remover"**



**FINE Red RDA 120**

		Usually	NOW
Tube 60ml	DA-7311-12	\$48.00	<b>\$40.80</b>
Paraben Free	DA-7311-15	\$48.00	<b>\$40.80</b>
Single Dose	DA-6901-21	\$215.00	<b>\$182.70</b>

**"Ideal for final polishing"**



**EXTRA FINE Yellow RDA 40**

		Usually	NOW
Tube 60ml	DA-7311-11	\$48.00	<b>\$40.80</b>
Single Dose	DA-6901-20	\$215.00	<b>\$182.70</b>

**"Clean, polish and desensitize"**



**CLEAN & POLISH Lilac**

		Usually	NOW
Tube 60ml	DA-7411-11	\$56.00	<b>\$47.60</b>
Single Dose	DA-6901-25	\$270.00	<b>\$229.50</b>

**"Powerful relief from dentin hypersensitivity"**




**PROPHY PASTE HAp**

		Usually	NOW
Tube 60ml	DA-7311-17	\$48.00	<b>\$40.80</b>



# Signature Instruments

Comfort, Quality & Durability



Pre-Molar Orca Black

Diagnostic Pacific Blue

Universal Nautical Blue

Anterior Passionate Purple

Distal Caribbean Green

Implant Tangerine Orange

Mesial/Posterior Sunshine Yellow

Restorative Alabaster Grey

Tissue Grafting Bermuda Blue

## Montana Jack Sickle



PT-R138

Contra angle for the most universal adaptation of all instrument designs.

## Montana Jack Rigid



PT-R138R

Adapts like the original Montana Jack but with extra backbone for the removal of tenacious calculus, cement and overhangs.

## Montana Jack Pink



PT-R138BC

Contra angle for the most universal adaptation of all instrument designs.

## Queen of Hearts



PT-R144

This periodontal curette has long fine cutting edges for finishing with a completely smooth surface when root planning, bringing your patient one step closer to better overall health.

## Jack B. Nimble



PT-R130

Small, contra-angled scaler with a continuous curved blade. The blade is structured so increased surface of the cutting edge comes in direct contact with the tooth surface when scaling.

## Jack B. Nimble Slim



PT-R130S

The Jack B. Nimble® is a small, contra-angled scaler with a continuous curved blade. The blade is slimmer than the PTR130.

## H5-L5



PT-R090

A '2 in 1' anterior instrument with an H5 straight sickle on one end for deposit removal and interproximal access, and a Langer 5 on the opposing end for fine finishing.

## H5-L5 Mini



PT-R092

A scaler and curette combined into one instrument adapting like the original H5-L5 above but the Langer 5 mini curette increases subgingival adaptation for the narrower anterior anatomy.

## MT 4-5



PT-R120

Adapts like the 204s with the same angulation but it is more rigid and has an extended reach for removal of tenacious calculus.

## NEB 128-L5 Mini



PT-R093

This design is very popular in pedo practice with the narrow Nebraska 128 sickle on one end and the Langer 5 on the other end.

## Mini me



PT-R089

A full anterior universal blade on one end and a mini anterior universal blade on the other end.

## MT Explorer



PT-R170

Same Contra-angle as the Montana Jack. Durable spring tip provides convenient interproximal access. Very thin sickle curvature.

## Mirror Handle



PT-R245

**Cone Socket.** Ideally weighted, large diameter, ergonomic cone socket mirror handle with mobility probe.

## Probe WRS ACE



PT-R215

**UNC12.** Advanced Comprehensive Exam (ACE™) Probe with 3-6-9-12mm markings and a straight probe end with 1-12mm markings.

## ACE Furcation Probe



PT-R216

**PR12 Short.** Advanced Comprehensive Exam Probe 3-6-9-12 mm markings on each end, one end standard probe, one furcation probe.

**Buy 4 Get 2 Free!**

All signature instruments Gracey, Scalers, Composite instruments, universal curettes. Excludes Wingrove, probes and explorers.

**\$87.00**

### Woodpecker Ultrasonic Scaler Handpiece EMS Fitting

Ultrasonic Scaler Piezo  
Handpiece Compatible with EMS

- High performance
- CE certified
- Autoclavable



**\$288.27**

Usually \$320.30

SC-OPAM-EMS

EMS  
LED Fitting

SCORPION  
Dental Instruments



SC-HW-5L

**\$288.27**

Usually \$320.30

10%  
OFF

**Buy 1 Scorpion Tip & Get 1 Wrench Free!**



\*Not compatible with new generation Kavo handpiece

CLiP

**IMPLANT SCALING TIPS**

SC-OI-EM

**\$161.20**  
ea

Scorpion tip adaptable on  
**EMS**



Scorpion tip adaptable on  
**Satelec /  
NSK piezo**



Scorpion tip adaptable on  
**Sirona  
Sirosonic**



Scorpion tip adaptable on  
**Kavo PiezoLux -  
SonoSoft Mectron**



Scorpion tip adaptable on  
**Kavo SONICflex  
2000-2003**



Scorpion tip adaptable on  
**W&H Synea-Proxeo,  
NSK Ti-Max S970,  
NSK AS2000**



**SCALING TIPS**

SC-OBEM

**\$139.70**  
ea

**PERIO SCALING TIPS**

SC-OPEM

**\$182.60**  
ea



# BEYOND® II

## Ultra Whitening Accelerator

The sleek modern design of the BEYOND® II Ultra Whitening Accelerator works for any dental office

The BEYOND® II Ultra whitening accelerator is a more portable option for in-office teeth whitening. The newest incarnation still allows for the simultaneous treatment of two patients. The tailored beam of blue light delivered by LEDs is only enhanced by the ultrasonic power that assists in providing superior whitening results in a short treatment time.

Variable settings ensure that each patient receives a comfortable, positive whitening experience specifically catered to them. The three output settings of high, medium, or low on the digital display makes this specialised treatment possible. The lightweight portable design enables travel and ease of use from room to room. Utilising negative ions, the built-in air purifier reduces airborne pollutants, actively contributing to a cleaner healthier environment.

**1**  
Remaining  
on hand



### Simultaneous

Allows for the effective treatment of two patients at the same time.

### Adjustable

Whitening option can be set to gentle, normal, or intense light output settings for a tailored treatment that best suits each individual patient.

### Digital

Simple and quick patient preparation and system operation with digital display control board.

### Elegant

Chic and modern design that will be easily embraced by any dental practice.

### Purifies

Using negative ions, the built in air purifier reduces airborne pollutants that may exist within the practice, actively contributing to a cleaner and healthier environment.

### Unique

Advanced design allows for a portable system, allowing for ease of mobility through different rooms within the practice.

**Includes 1 x BEYOND Five-Patient  
Economy Kit (BY-PD205) !!!**  
Valued at \$480.00

BY-BE-818UL

BEYOND II ULTRA WHITENING ACCELERATOR

**\$2,822.00**  
ea



## BEYOND® COMPLETE®

### Professional Teeth Whitening Kit with Tray Whitening Material

#### Single Patient, Comprehensive Whitening Package



This 35% hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) In-office teeth whitening kit includes everything needed for one patient.

**Kit Contents:** 1 Solo gel Dual-barrel syringe 35% (4g/3.5 mL), 2 syringes CoreWhite gel (mint) 6% (2.36g/2.2 mL), 1 BlueSeal Single syringe (2g/1.6mL), 1 Mixing tip, 1 Plastic tip, 1 Metal tip, 1 Cheek Retractor, 1 Jar of lip cream (2.4g/2mL), 1 Suction tip, 2 Sheets EVA impression tray material, 1 Face Protector, 1 Complete instruction manual, 1 Vita Shade Guide, 1 After-treatment guide

BY-PD202

Normally \$188.00

**\$169.20** ea

## BEYOND® SOLO®

### Single Treatment Kit

#### Our Most Economical, Single Patient Whitening Kit

This 35% hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) In-office teeth whitening kit includes the essentials for chairside teeth bleaching of one patient.



**Kit Contents:** 1 Solo gel dual-barrel syringe (2.6g/2.2 mL), 1 BlueSeal single syringe (2g/1.6mL), 1 Mixing tip, 1 Plastic tip, 1 Metal tip, 1 Face Protection Cloth, 1 Solo instruction manual, 1 Vita Shade Guide, 1 After-treatment guide

BY-SL101

Normally \$83.00

**\$74.70** ea

## BEYOND® Five-Patient Economy Kit

#### Our Most Economical, Full-sized Whitening Kit

This 35% Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>), In-office teeth whitening kit includes material for five patients.



**Kit Contents:** 5 Beyond single dual-barrel syringes (4g/3.5 mL), 5 BlueSeal single syringe (2g/1.6mL), 5 Mixing tips, 5

Plastic tips, 5 Metal tips, 5 Cheek Retractors, 5 Jars of lip cream (2.4g/2mL), 5 Face Protection Cloth, 1 Five-Patient instruction manual, 5 Vita Shade guides, 5 After-treatment guides

BY-PD205

Normally \$480.00

**\$432.00** ea

## BEYOND® OSMO®

### Five-Patient Professional Teeth Whitening Kit

#### 6% H<sub>2</sub>O<sub>2</sub> Premium Whitening Formula with Post-treatment Whitening Pen



This 6% Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>) Teeth Whitening Kit is designed to deliver professional results in just 36 minutes.

**Kit Contents:** 5 Osmo gel dual-barrel syringe (2.6 g / 2.2 ml), 5 Osmo teeth whitening pens (3.2 g each), 5 Mixing tip, 5 Face Protection Cloth, 1 Osmo instruction manual, 5 Vita Shade Guide, 5 After-treatment guide

BY-OS101

Normally \$376.00

**\$338.40** ea




# Toothbrushes

TePe's user-friendly toothbrushes are developed for efficient and gentle plaque removal. A variety of models and sizes facilitate daily oral care for adults and children.

We care for healthy smiles

## TePe Nova



- Easy-reach tip for access to hard-to-reach areas, especially around the back teeth
- Round-ended filaments ensure gentle cleaning
- Easy-grip handle with thumb pad for extra comfort



**\$51.00** ONLY \$3.64 per Brush (Blister Pack 14)  
**\$230.50** ONLY \$2.88 per Brush (Flow Pack 80)

## TePe Supreme



- Efficient cleaning due to the longer, tapered bristles
- Unique two-level filament design for excellent reach
- Easy-grip handle with thumb pad for extra comfort

**\$52.00** ONLY \$3.71 per Brush (Blister Pack 14, TP-352688)  
**\$258.10** ONLY \$3.22 per Brush (Flow Pack 80, TP-352100)

## TePe Supreme Compact



- Unique two-level filament design for excellent reach
- Easy-grip handle with thumb pad for extra comfort
- Suitable for children or for adults who prefer a smaller brush head

**\$258.10** ONLY \$3.22 per Brush (Flow Pack 80, TP-392120)

## TePe Special Care



- First brush stage after oral surgery
- For patients with extremely sore and delicate tissues, as a result of trauma and radiotherapy treatment

**\$61.20** ONLY \$4.37 per Brush (Blister Pack 14, TP-452357)  
**\$95.90** ONLY \$3.83 per Brush (Flow Pack 25, TP-452150)

## TePe Special Care Compact



- First brush stage after oral surgery
- For patients with extremely sore and delicate tissues, as a result of trauma and radiotherapy treatment

**\$95.90** ONLY \$3.83 per Brush (Flow Pack 25, TP-452151)

## TePe Implant Orthodontic



- Long and thin neck for better access
- For cleaning buccal implant surfaces or orthodontic appliances

**\$54.10** ONLY \$3.86 per Brush (Blister Pack 14, TP-442368)  
**\$75.50** ONLY \$3.02 per Brush (Flow Pack 25, TP-442100)

## TePe Extra Grip



- Provides a comfortable, stable grip which helps patients with reduced manual dexterity
- Weighs only 30 grams
- Fits TePe Select, TePe GOOD toothbrushes and most of our special brushes

**\$40.80** ONLY \$10.20 per Grip (Pack 4, TP-912210)

# Specialty and Kids Toothbrushes

TePe offers a wide range of cleaning devices developed for individuals with specific oral hygiene needs.



We care for healthy smiles

## TePe Compact Tuft



- Dome-shaped brush head with dense rounded tuft
- Ideal for cleaning along the gumline, on the occlusal surfaces of emerging teeth, or impacted wisdom teeth



TP-432368

**\$62.20**

ONLY \$4.44 per Brush



TP-432100

**\$82.60**

ONLY \$3.30 per Brush

## TePe Universal Care



- A unique angle provides easy access to the lingual and palatal surfaces of teeth and implants
- Facilitates cleaning in narrow and hard-to-reach areas
- End-rounded filaments to ensure a safe and gentle clean



TP-472368

**\$65.30**

ONLY \$4.66 per Brush

## TePe Interspace



Soft

TP-412364



Medium

TP-412366

- Cleaning behind furcations, implants, fixed braces or the last tooth
- Brush tip can be attached to either side of the neck for different angles
- Comes with 12 replaceable brush tips
- Colour-coded tip softness choices

Yellow Tips: soft  
Blue Tips: medium



**\$118.30**

ONLY \$8.45 per Tip

## TePe Kids



- Brush with colourful print suitable from the age of 3 years
- Tapered small brush head for easy access
- Round-ended filaments for gentle and efficient cleaning
- Handle suits both kids and assisting adults



Select Compact with illustrations.

TePe Kids are recommended from 3 years and up.

Available in Extra Soft or Soft.

Extra Soft TP-339160  
Soft TP-339161



**\$83.60**

ONLY \$3.34 per Brush

## TePe Mini



- Recommended from the first tooth
- Tapered small brush head and round-ended filaments for efficient cleaning
- Handle provides a stable grip for the adult brushing the child's teeth, but fits just as well in a small hands



Recommended from the first tooth.

Available in Extra Soft or Soft.

Extra Soft TP-382110  
Soft TP-382120



**\$51.00**

ONLY \$2.04 per Brush

# GERMIPHENE®



## Ultra Swipes Wipes Refill

### Hard Surface Disinfectant Towelettes

Ultra Swipes Plus ensures broad spectrum disinfection in both the presence and absence of whole blood against Poliovirus Type 1 and TB (Mycobacterium bovis).

160 Pack **\$24.00**  
GP-HSWR

Buy 12 or more **10% OFF**  
Buy 24 or more **20% OFF**



## Gamacide 3 Wipes Refill

### Multi-Surface Disinfectant Cleaner

Gamacide3 is the ideal surface disinfectant for your cleaning and disinfecting needs. Gamacide3 kills bacteria, non-enveloped difficult to kill viruses such as polio, TB and fungi.

160 Pack **\$24.00**  
GP-HGAM3W

Buy 12 or more **10% OFF**  
Buy 24 or more **20% OFF**



## Genie Plus Wipes Refill

### Chair Disinfectant / Cleaner

ALCOHOL FREE. Cleans and polishes naugahyde, plastic laminate, vinyl and leather with no sticky or slippery residue. It removes built-up dirt other solutions leave behind.

160 Pack **\$22.00**  
GP-HGENWIPE

Buy 12 or more **10% OFF**  
Buy 24 or more **20% OFF**



## Gamacide 3 Solution

### Multi-Surface Disinfectant Cleaner

Gamacide3 is the ideal surface disinfectant for your cleaning and disinfecting needs. Gamacide3 kills bacteria, non-enveloped difficult to kill viruses such as polio, TB and fungi.

4L Bottle **\$88.20**  
GP-HGAMA3  
Usually \$98.00



## Germxtra Hard Surface Disinfectant

- Rapid, effective ready-to-use disinfectant spray
- Bactericidal, Fungicidal, Virucidal and Tuberculocidal
- Neutral pH formula
- Synergistic blend of ethanol (ethyl alcohol) and surfactants that clean and disinfect in one step

3.78L Bottle **\$111.00**  
GP-HGERXX4  
Usually \$123.40



## Gleam II Ultrasonic Cleaning Solution

- Makes your instruments shine
- Will not harm stainless steel instruments
- Safe on hard plastics, glass and rubber
- Easy to use 2L Tip N' Measure bottle
- 2L Gleam II = 22L solution

2L Bottle **\$100.90**  
GP-HGLE  
Usually \$112.20



## Gzyme3 Enzymatic Cleaning Solution

- Penetrates and removes organic matter from all surgical and dental instruments
- Quickly loosens and dissolves proteinaceous deposits
- Bacteriostatic to prevent additional microbial growth
- Neutral pH formula is mild on instruments
- Non-corrosive, low foaming, free-rinsing
- 1L Gzyme3 = 400L solution

## More for less...

- 1Litre of Gzyme3 Makes 400 Litres of Solution
- Super concentrated, Super economical
- Uses little storage space
- 1 Bottle to recycle, versus 400 bottles
- Convenient Squeeze & Measure packaging = Less mess, no spills, no waste



Tangerine Scented  
1L Bottle **\$155.10**  
GP-HGZY3  
Usually \$172.40



**E20**  
Dental Waterline  
& Shock Treatment



An easy and effective way to control biofilm and bacteria in your water system, without using harsh chemicals. Reduce the risk of cross contamination, while cleaning your dental waterline tubing. Specifically designed for use in independent water bottle systems.

236ml  
Bottle  
GP-HE20

**\$119.20**  
Usually \$149.00



**Microvac**  
Enzymatic Evacuation  
System Cleaner

Microvac is an effective multipurpose cleaner with a unique low foaming multi-enzyme formula to effectively remove organic debris, bioburden and proteinaceous materials. Use Microvac to safely clean rubbers, plastics, and stainless steel.

4 Litre Bottle  
(no pump)  
GP-HMICRO

**\$77.10**  
Usually \$85.70



**Gluterate**  
Sterilising and  
Disinfecting Solution

- 2% Glutaraldehyde solution
- Ready to use, requires no activation
- Bactericidal, Fungicidal, Tuberculocidal, Virucidal, Sporocidal
- 10 Hour sterilisation
- 10 Minute high level disinfection

3.78L  
Bottle  
GP-HGLUT

**\$84.40**  
Usually \$93.80



**Glute Neut**  
Glutaraldehyde Neutralizer

- For the safe disposal of Glutaraldehyde
- Fast-acting neutralizer
- Eliminates vapours, eye and skin irritation
- Changes solution from yellow to red to assure that neutralization has occurred

16 x 50g  
Packets  
GP-HGLUTNEUT

**\$131.30**  
Usually \$145.90



**IC Foam**  
Fragrance Free  
Hand Sanitiser Foam

- One step sanitising and moisturising
- 70% Ethyl alcohol
- Tough on germs, gentle for repeated use
- Stable, fluffy foam that applies like cream
- Fast drying for glove donning
- Non sticky, non greasy
- Dye free, non staining, non aerosol foam

550ml  
Pump Bottle  
GP-HICFOAM

**\$27.50**  
Usually \$30.60



**IC Gel**  
Antiseptic Skin Gel

- Broad spectrum
- 70% Ethanol + emollients to moisturise skin, made with chamomile
- Designed for use between patient contact to reduce cross-contamination
- Dye free, non staining, residue free
- Evaporates quickly to facilitate glove donning
- Does not impact glove performance

2Litre  
Bottle  
GP-HIC2L

**\$60.50**  
Usually \$67.30

480ml  
Pump Bottle  
GP-HIC

**\$22.90**  
Usually \$25.50



**E.T. Foam**  
Enzymatic Transportation Foam

- Ready to use, instrument pre-cleaning solution
- Breaks down bioburden to speed up cleaning process
- Keeps soiled instruments moist
- Ideal for use prior to instrument washers – rinse prior to placing in instrument washer
- Neutral pH, non-corrosive formula will not harm most metals, plastics, rubber or corrugated tubing



1L Spray  
Bottle  
GP-HET

**\$55.20**  
Usually \$69.00

Virofex®  
**ViroVent**<sup>365</sup>

**Aspirator Disinfectant**  
**Daily Concentrate**

Highly effective and economical disinfectant concentrate used at only 2% concentration. Bactericidal, Yeasticidal, Fungicidal, Tuberculocidal and effective against viruses including: HBV, HCV, and HIV.

1L  
 TD-VIR0018  
**\$41.30**  
 Usually \$45.90



Virofex®  
**ViroSurf**

**Surface Disinfectant**  
**Alcohol-Based**

Highly effective surface disinfectant designed for use on non invasive medical device surfaces. Bactericidal, Yeasticidal, Fungicidal, Tuberculocidal, Mycobactericidal and effective against viruses.


**\$82.60**  
 Usually \$91.80 5L  
 TD-VIR0032



Virofex®  
**ViroSurf**  
 wipes

**Surface Disinfectant Wipes**  
**Alcohol-Based**

Highly effective surface disinfectant wipes designed for use on non-invasive medical device surfaces. Bactericidal, Yeasticidal, Fungicidal, Tuberculocidal, Mycobactericidal and effective against viruses including Coronavirus, HBV, HCV, HIV, Adenovirus and Norovirus.



Buy  
**12**  
 or more

**25%**  
 OFF

Refills  
 160 Wipes  
 TD-VIR0028  
**\$16.50**  
 Usually \$22.00

Empty Tub  
 TD-VIR0028  
**\$7.70**

### Cotton Rolls

High quality, soft, made from 100% virgin cotton, are non-linting, non-sterile, pliable and 20% more absorbent than the leading brand.



DF-CS-0200 2000 /Box

Normally \$33.70 **\$30.30** ea

### Non-Woven Gauze

Made from a 4-ply, 32 gram rayon/polyester blend which provides for maximum absorbency and little adhesion to wounds.



DF-NW-0200 2" x 2" x 200

Normally \$6.10 **\$5.40** ea

On Backorder. Stock arriving early 2024

### Dry Dent

Soft, comfortable and highly efficient

### DIRECTA

**\$44.00** Small **\$56.00** Large



**\$26.00**

Sublingual

Small Parotid x50	DA-4160-73
Large Parotid x50	DA-4160-72
Small Sublingual x50	DA-4160-69
Large Sublingual x40	DA-4160-68

### Dry Tips



Small x50	MA-291543
Large x50	MA-291542
Small Reflective	MA-291727
Large Reflective	MA-291726

**\$28.60** Regular **\$39.80** Reflective Small **\$43.90** Reflective Large

## Microbrush International

### Tube Series 400 Packs



**\$55.10**

Regular 2.0 mm	MA-400
Fine 1.5 mm	MA-F400
Superfine 1.0 mm	MA-SF400

### Regular 400 Packs



**\$55.10**

Blue	2.0 mm	MA-MRB400
Green	2.0 mm	MA-MRG400
Purple	2.0 mm	MA-MRP400
Peach	2.0 mm	MA-MRS400

## The Original Microbrush® Applicator

Consistent, reliable quality

### Plus Series 400 Packs



**\$55.10**

Super Fine-White	MA-PSF400
Fine Pink	MA-PF400
Regular Blue	MA-PR400

### Extended Series 100 Packs

Extended Reach Refill MA-PX100



**\$35.70**

Extended Reach with dispenser MA-PX50




**\$34.70**

### Surgical Evacuation Tubes 6"

For micro-aspiration use White (1/8")  
For macro suction use Green (1/4") or Blue (3/8")

- Pre-bent to a 30° angle for greater access
- AUTOCLAVABLE



White 1/8in  
25 Pack  
V-403201

Green 1/4in  
25 Pack  
V-403202


Blue 3/8in  
25 Pack  
V-403211

**\$35.00**

**\$28.80**

Usually \$39.00

Usually \$32.00



### Surgical Aspirator Tips

Made with smooth edges for maximum patient comfort. The surgical aspirator is molded at a 30 degree angle and fits 11mm suction holders.

Blue 1/16in  
25 Pack  
DF-ST-1020

White 1/8in  
25 Pack  
DF-ST-1021

Green 1/4in  
25 Pack  
DF-ST-1023

**\$10.90**

Usually \$12.20

### Saliva Ejector 15cm

- Transparent with fixed tip
- Fully flexible to bend to the required angle
- Phthalate Free



Blue  
100 Pack  
CT-9008272

**\$8.90**

Usually \$9.90



### Hygofix Saliva Ejectors

- Prebent saliva ejector with integrated tip to fit 6.5 mm system
- Manufactured in one piece and one material
- Great if you don't have an assistant, stays in place
- Colour: White

100 Pack  
OR-HF1100

**\$23.80**

Usually \$26.50

### Patient Bibs

2-ply tissue and 1-ply poly for absorption and patient protection. 13" x 18"

Blue  
Case 500  
DF-PB-8001

**\$45.90**

Usually \$51.00



### Patient Bibs

2 layers of absorbant paper tissue and 1 layer of polyfilm. 33 x 45cm.

Blue  
Case 500  
AD1-ADA201

**\$35.80**

Usually \$39.80



### Dappen Dishes



Assorted  
1000 Pack  
PP-004B

**\$80.60**

### Drinking Cups




Blue  
1000 Pack  
DF-DC-7001

**\$50.40**

Usually \$56.10

### Micro-Stix

- Smart hold technology lets the practitioner determine the holding strength
- Bendable neck for hard-to-reach areas
- No adhesive residue left on surfaces



64pk Original  
Hold Blue  
MA-STIX64B

**\$48.60**

Usually \$54.10

### Pic-n-Stic

The original "handle for small objects" has numerous applications from dentistry to model making to replacing hearing aid batteries.



64pk Original  
Hold Blue  
PL-DPIC

**\$33.30**

Usually \$40.00



### Headrest Covers

These Headrest Covers are designed to fit over the top of headrests. They are made of clear plastic and are easy to apply and dispose of.

9.5" X 11" ( 24 x 28cm)  
Clear Box of 250

DF-BF-9000

9.5" X 14" ( 24 x 35.5cm)  
Clear Box of 250

DF-BF-9500

**\$29.30**

Usually \$32.60



### Light Handle Sleeves

These Light Handle Sleeves are made of clear plastic and are easy to apply. Made to fit most light handles.

4" x 5.75"  
Clear Box of 500

DF-BF-5000

**\$21.10**

Usually \$23.50



### Syringe Sleeves Air/Water



2.5' x 10" (6.3 x 25cm)  
Clear Box of 50

DF-BF-3000

**\$19.20**

Usually \$21.40

### Syringe Barrier Sleeves

Protective Disposable Barrier Sleeves created to control cross infection for all flowable and injectable syringes.



Clear Box of 50

SH-2070

**\$7.30**

Usually \$8.20



### Curing Light Sleeves

FlashMax P3 Disposable Sleeves.

For use with the FlashMax P3 Curing Light. Disposable and a hygienic way to avoid cross contamination.



Clear Box of 250

CM-S100513

**\$33.90**

Usually \$37.70

### Curing Light Sleeves

Made from very supple poly plastic for greater tactile feel. Easy to apply.



Pen Type Small  
Clear Box of 500

PP-123S

Pen Type XS  
Clear Box of 500

PP-123SS

**\$43.10**

Usually \$47.90

### Barrier Film

- Perforated 4" x 6" sheets leave behind no residue.



Barrier Film Clear NS Edge 1200 Boxed

DF-BF-1500

Barrier Film Blue NS Edge 1200 Boxed

DF-BF-2500

**\$38.50**

Usually \$42.80

Barrier Film NS Edge  
No Box 1200

DF-BF-2600

**\$33.30**

Usually \$37.00

### Flecta Mirrors

- A shiny new mirror for every patient
- The highest quality disposable mirror available
- 40% Larger viewing surface

**PULPDENT®**



PL-DFLEC72

**\$93.60**

Usually \$104.00

### Cyber Tech Level 2 Earloop

Softer inner lining of the mask provides greater comfort and allows for better breathability

Box of 50  
CT-9002733

Blue  
CT-9002783  
Pink

**\$8.10**

Usually \$9.00



### PriMED 1273 Level 3

- Highest level of fluid protection based on ASTM F2100-11 standards
- Latex free
- Anti fog

Box of 50  
PS-SZ-PG4-1073

**\$22.95**

Usually \$25.50



### PRIMAGARD 306 Level 2

- Fluid resistant to 16kPa (120mmHg)

Box of 50  
PS-PM4-306

**\$12.78**

Usually \$14.20



### Level 3 Dual Fit Ear-Loop Face Masks (Pleated)

- Fluid resistant outer layer
- 4 Ply
- Dual fit chin contour technology
- Soft, rounded ear loops

Box of 50  
DF-MK7300-C

**\$22.50**

Usually \$25.00



### Mediflex Level 2 Earloop

- Wand (Medsafe NZ) registered
- Pleated
- Form fit

Box of 50  
EZ-MUP102

**\$18.20**

Usually \$20.30



### Full Face Visors & Eye Shields

- Lightweight & comfortable to wear
- Optical grade lens

The anti-fog, latex free visors are designed to protect your eyes and face by preventing exposure from hazardous or infectious substances.



Description	Code	Usually	Now Only
Full Face Visor Starter Pack 3 Visors and 15 shields	MS-VMED1000-03	\$144.80	<b>\$130.30</b>
Full Face Visor Shield	MS-VMED1000-06	\$183.60	<b>\$165.20</b>
Eye Shield Frames Assorted 10 Pack	MS-VMED100-80	\$59.20	<b>\$53.20</b>
Eye Shield Lenses 50pk	MS-VMED100-92	\$203.00	<b>\$182.70</b>



### Nitrile Powder-Free Examination Gloves Blue

Extra Small	CT-9889873
Small	CT-9889874
Medium	CT-9889875
Large	CT-9889876

Box Normally **\$8.90**  
100 \$9.90



### Latex Powder-Free Examination Gloves Natural latex colour

Extra Small	CT-9884705
Small	CT-9884706
Medium	CT-9884707
Large	CT-9884708

Box Normally **\$13.50**  
100 \$15.00



### Luer Lock Syringes



**10% OFF**

#### Dripless Luer Lock Syringe 6cc Irrigation Tips 100 Pack

V-318060



**\$106.20**

Usually \$118.00



Luer Lock Syringe  
100 Packs

	Usually	10% OFF
1.2cc V-316002	\$91.00	<b>\$81.90</b>
3cc V-316030	\$42.00	<b>\$37.80</b>
6cc V-316060	\$58.00	<b>\$52.20</b>
12cc V-316012	\$58.00	<b>\$52.20</b>

#### Appli-Vac Pre-Tipped 3cc Luer Lock Syringe 100 Packs



23G V-316133	<b>\$90.00</b>
27G V-316137	
30G V-316139	

Usually \$100.00



## SWEEP AG200 ULTRASONIC CLEANER With Basket

The Sweep 200 with SweepZone® Technology delivers uniform ultrasonic cleaning throughout its 1.25 gallon (4L) tank.

### Designed and Manufactured in the USA.

- Digital touchpad timer with 5 preset cleaning cycles for easy operation and automatic shutoff
- SweepZone Technology delivers powerful and repeatable cleaning results
- Multi-positional drain
- Designed for countertop use

### Specifications

- Input Power 117V 60Hz: 100 Watts
- Amps: 0.8 Amps
- Output/Cleaning Frequency: 43kHz (+/- 2kHz)
- Tank Capacity: 4.0 L (1.06 Gal)
- Internal Tank Dimensions: 29.8 x 15.2 x 10.2 cm
- Overall Dimensions: 32.2 x 17.5 x 22.2 cm

LR-AG200


Usually \$2,340.90

**\$2,106.<sup>80</sup>**

Constructed of stainless steel, the Sweep 200 features a digital touchpad timer with 5 preset cleaning cycles ranging from 6 to 60 minutes. Operating at multi cleaning frequencies, 43kHz +/- 2 kHz, SweepZone Technology eliminates hot spots and weak areas throughout the ultrasonic cleaner by uniformly sweeping ultrasonic power across the tank. Using Pulse width Modulation (PWM), the Sweep 200 senses variations in cleaning loads, solution levels, and voltage fluctuations to deliver consistent and powerful cleaning performance.



## WAVE CHECK TEST



WaveCheck® is a specially designed ready-to-use indicator designed to test the efficacy of your ultrasonic cleaner. Now more than ever, instrument reprocessing is a critical part of maintaining patient and staff safety. Instruments need to be thoroughly cleaned and decontaminated prior to sterilisation.

During the test, the monitor strip is measuring multiple variables of the ultrasonic cleaning process, including cavitation, time, temperature, solution ratio and ability. Full removal of the test soil from the strip indicates an effective cleaning process.

LR-WAVECHECK

**50 Test Strips Per Pack**  
**Free Test Strip Holder included**


Usually \$255.00

**\$229.<sup>50</sup>**



## Sterilisation Pouches with Built-in Internal Indicator

- Thumb notch for easy opening of pouch
- Heavily sealed and reinforced bottom for added protection against instrument breaching
- Triple heat-sealed pouch to prevent possible penetration and tearing from sharp instruments
  - Medical grade paper that meets or exceeds medical pouch standards
  - Built in internal and external indicators
- Transparent film pouch for easy viewing of contents



**Pre folded flap helps provide a tight seal!**



DEFEND.LOC sterilisation pouches are manufactured with lead-free, built-in dual internal and external indicators to ensure the correct sterilisation temperature is reached in the autoclave chamber as well as inside the instrument compartment. Eliminates the need for internal indicator strips.

**DF-SP-0150** 2.25" X 4" (5.7 x 10.16 cm) 200 /BOX

**DF-SP-0550** 2 3/4" X 10" (7 x 25.4 cm) 200 /BOX

**DF-SP-1500** 3 1/2" X 5 1/4" (9 x 13.4 cm) 200 /BOX

**DF-SP-2500** 3 1/2" X 10" (9 x 25.4cm) 200 /BOX

Normally \$20.40

**\$18.30** ea

# Happy Holidays!

As the sun sets on another year, we would like to thank you for your ongoing business and support throughout 2023.

We look forward to working with you again throughout 2024.

*Wishing you a safe, fun and relaxing holiday!*

*The Team at*



**Closing 12pm Friday 22<sup>nd</sup> December 2023 // Reopening Wednesday 3<sup>rd</sup> January 2024!**